

supporty[®]
Orthopedic Line

PRODUCT CATALOGUE



Partner[®]

Let Us Improve Your Health

PARTNER, Mobility Aids, orthopedics designer & manufacturer since 2014.

Our Vision is to be a major player in the world of medical equipment. The PARTNER has specialized over the years in the manufacturing and designing of orthopedics, wheelchairs, Commode chairs, comfort chairs, Crutches, Canes, Rollator and any other product that may accompany a patient or user in maintaining its autonomy and mobility.

Our European development and design teams are equipped with precision high-tech tools resulting in reliability and quality for each product we produce.

Today, PARTNER is a group of more than 100 employees with subsidiaries in Europe and Asia and delivers products sold in more than 15 countries worldwide.



C O N T E N T S

CERVICAL COLLARS

01-03

SHOULDER SUPPORTS

04-07

ELBOW SUPPORTS

08-09

WRIST AND FINGER SPLINTS

10-17

CORSETS

18-23

KNEE SUPPORTS

24-31

ANKLE SUPPORTS

32-36

FOOT CARE

37-45

PEDIATRIC PRODUCTS

46-52

OTHER PRODUCTS

53-55

CERVICAL COLLARS



CERVICAL COLLARS

CERVICAL COLLAR (Sponge)

Made of compressed Sponge witch is covered by cotton fabric.
Keeps the neck warm and stabilize the neck and the head.

USAGE

Muscles pain.
Neck arthritis.
Sliding discs.
After the traumas and bruises



SU-1001

Size:

S 34-36 cm
M 36-38 cm
L 38-40 cm
XL 40-42 cm



NELSON COLLAR



Made of the foam witch has holes in it that helps coaling the neck, allows the x-ray to goes throw it, non-allergic, reduces the pain by reducing the ability of the head and the neck to move.

USAGE

Muscles pain.
Neck arthritis.
Sliding discs.
After the traumas and bruises

SU-1002

Size:

S 34-36 cm
M 36-38 cm
L 38-40 cm
XL 40-42 cm



PHILEDELPHIA COLLAR

Made of the foam, reduces the ability of moving the neck and head thanks to the plastic supports on the front or the back, allows the x-ray to go throw it.

USAGE

Soft tissue trauma.
Ligaments and bones injuries.
Sliding discs.
Stabilizations in the cases of the uncon-
sciousness.



SU-1003

Size:

S 34-36 cm
M 36-38 cm
L 38-40 cm



PHILEDELPHIA COLLAR (Front Open)



Made of the foam, reduces the ability of moving the neck and head thanks to the plastic supports on the front or the back, allows the x-ray to go throw it, the front has a hole in the Trachea area to be used in case of emergency.

USAGE

Soft tissue trauma.
Ligaments and bones injuries.
Sliding discs.
Stabilizations in the cases of the uncon-
sciousness.

SU-1004

Size:

S 34-36 cm
M 36-38 cm
L 38-40 cm
XL 40-42 cm



CERVICAL COLLAR WITH / WITHOUT CHIN SUPPORT (Plastic)

Made of PVC with chin support in a design that can be modified as wanted to stabilize the neck and reduce the pressure on the spine

USAGE

Muscles pain.
Neck arthritis.
Sliding discs.
Cases of radiculopathy
After the traumas and bruises



SU-1005 / SU-1006

Size:

S 34-36 cm
M 36-38 cm
L 38-40 cm
XL 40-42 cm



WIRE-FRAME COLLAR JR



Made of the aluminum covered by polyurethane to stabilize the neck as wished.

USAGE

Discs hernia.
Arthritis.
Recovery from the surgery.
Neck injuries.
Tissue injuries

SU-1007

Size:

S 34-36 cm
M 36-38 cm
L 38-40 cm
XL 40-42 cm



EMERGENCY NECK COLLAR

Made of the PVC plastic with a foam layer on the inner face of it.
Can be adjusted to 3 different sizes.
Has a lock system that helps stabilizing the neck.
Has a hole in front of it to provide faster access to the neck in cases of emergency.
Allows the X-ray to go through it.

USAGE

Stabilizing the neck during the emergencies.
In coma cases.



SU-1008

Size:

S 34-36 cm
M 36-38 cm



MINERVA ORTHOSIS



Made of the Thermoplastic with is covered from the inside by sponge that help stabilizing the neck and the upper part of the body.
Stabilize the chest and reduces the neck's ability of movement.

USAGE

After the surgeries of the muscles and ligaments.
Arthritis.
Bone fractures.

SU-1009

Size:

S 34-36 cm
M 36-38 cm
L 38-40 cm



SHOULDER SUPPORTS



ARM SLING (Mesh Fabric Black)

Made of the 100% Polyester, nice on the skin, two layers of fabric provides more stabilization for the shoulder and the elbow and reduces the pressure on them.

USAGE

Shoulder dislocations
After surgeries
prevent edemas
Fractures



SU-2001

Size:

S 34-38 cm
M 38-44 cm
L 44-50cm
XL 50-58 cm
2XL58-66cm
3XL66-72cm



ARM SLING LUX



Made of sponge that is covered with the terry fabric from the inside. Stabilize the elbow and the shoulder and reduces the pressure on them.


USAGE

Shoulder dislocations
After surgeries
prevent edemas
Fractures

SU-2002

Size:

S 34-38 cm
M 38-44 cm
L 44-50cm
XL 50-58 cm
2XL58-66cm
3XL66-72cm



ARM SLING GABARDIAN

Made of Gabardian fabric. Stabilize the elbow and the shoulder and reduces the pressure on them.

USAGE


Fractures and bruises of the arm.
Bruised shoulder.



SU-2003

Size:

S 34-38 cm
M 38-44 cm
L 44-50cm
XL 50-58 cm
2XL58-66cm
3XL66-72cm



VELPAU BANDAGE UNIVERSAL



Made of sponge that is covered with the terry fabric from the inside. Stabilize the shoulder and the upper side of the arm. The straps around the chest provides stabilization in case of dislocations

USAGE

Shoulder dislocations
After arm and shoulder surgeries
prevent edemas
Fractures

SU-2004

Size:

UNISIZE

PADDED ARM SLING-(45°-60°)

Made of sponge that is covered with the terry fabric from the inside.
Stabilize the shoulder and the upper side of the arm on 45°-60°.
The straps around the chest provides stabilization in case of dislocations

USAGE

Shoulder dislocations
After surgeries
prevent edemas
Fractures
Armpits infections



SU-2005

Size:

UNISIZE

PADDED ARM SLING-30°



Made of sponge that is covered with the terry fabric from the inside.
Stabilize the shoulder and the upper side of the arm on 30°.
The straps around the chest provides stabilization in case of dislocations

USAGE

Shoulder dislocations
After surgeries
prevent edemas
Fractures

SU-2006

Size:

UNISIZE

CLAVICULA BANDAGE

Made of sponge that is covered with the terry fabric from the inside.
Stabilize the clavicle thanks to the straps that pulls the shoulder to the back to give the appropriate position.

USAGE


Clavicle dislocations
Clavicle fractures



SU-2007

Size:

S 70-80 cm
M 80-90 cm
L 90-100 cm
XL 100-110 cm
2XL 110-120cm



CLAVICULA BANDAGE



Made of sponge that is covered with the terry fabric from the inside.
Stabilize the clavicle thanks to the straps that pulls the shoulder to the back to give the appropriate position.

USAGE

Clavicle dislocations
Clavicle fractures

SU-2008

Size:

S 70-80 cm
M 80-90 cm
L 90-100 cm
XL 100-110 cm
2XL 110-120cm



CLAVICULA BANDAGE

Made of sponge that is covered with the terry fabric from the inside. Stabilize the clavicle thanks to the straps that pulls the shoulder to the back to give the appropriate position.

USAGE

Clavicle dislocations
Clavicle fractures



SU-2009

Size:

S 70-80 cm
M 80-90 cm
L 90-100 cm
XL 100-110 cm
2XL 110-120cm



SHOULDER SUPPORT (Neoprene)



Made of perforated neoprene fabric in 3 mm thickness. Provides warmth to the shoulder with protection in cases of the minor pain.

USAGE

Shoulder dislocations
After surgeries
Warming the Shoulder

SU-2010

Size:

S 34-38 cm
M 38-44 cm
L 44-50cm
XL 50-58 cm
2XL 58-66cm



SHOULDER SUPPORT (Wool)

Made of elastic cotton fabric specially designed to reduce sweating. Its inside is covered by smooth non-allergenic fabric. Designed to fit under clothes.

USAGE

Warming the Shoulder



SU-2011

Size:

S 34-38 cm
M 38-44 cm
L 44-50cm
XL 50-58 cm
2XL 58-66cm



SHOULDER ORTHOSIS



Made of thermoplastic covered by terry fabric from the inside. Has an optional elastic ligament. Stabilize the shoulder and the upper side of the arm.


USAGE

Stabilizing the shoulder and the upper part of the arm.
After surgeries.
Protection and Stabilization.

SU-2012

Size:

S 32-35 cm
M 36-39 cm
L 40-43 cm



ELBOW SUPPORTS



TENNIS ELBOW SUPPORT

Made of neoprene fabric that provides warmth and protection. Thanks to the three-headed silicone pad inside, it applies constant pressure and compression to the epicondyle heads.

USAGE

Lateral epicondylitis.
Internal articular fluid.
Elbow pain



SU-3001

Size:

S	22-25 cm
M	25-28 cm
L	28-31 cm
XL	31-34 cm
2XL	34-37 cm



EPICONDYLITIS BANDAGE UNISIZE

Made of neoprene fabric that provides warmth and protection. Thanks to the three-headed silicone pad inside, it applies constant pressure and compression to the epicondyle heads.

USAGE

Used in the treatment of epicondylitis by decreasing the load on the tendons by applying compression.



SU-3002

Size:

UNISIZE

PADDED ELBOW SUPPORT

Made of neoprene fabric that provides warmth and protection.
Used by athletes in cases of minor pains

USAGE

Protects elbow area against impacts.



SU-3003

Size:

S	22-25 cm
M	25-28 cm
L	28-31 cm
XL	31-34 cm
2XL	34-37cm



TENNIS ELBOW SUPPORT UNISIZE



Made of thermoplastic covered by terry fabric from the inside.

The angle can be adjusted easily
The Upper part can be adjusted in (0-15-30-45-60-75-90-105) degrees
The Lower part can be adjusted in (0-15-30-45-60-75-90-105) degrees

USAGE

Enables controlled mobilization of elbow.
medial and lateral epicondyle fractures.
After Surgeries.
elbow dislocation,elbow arthrolysis.

SU-3006

Size:

UNISIZE

ELBOW CONTRACTURE ORTHOSIS

Made of thermoplastic covered by terry fabric from the inside.

The angle can be adjusted easily
The Upper part can be adjusted in (0-15-30-45-60-75-90-105) degrees
The Lower part can be adjusted in (0-15-30-45-60-75-90-105) degrees

USAGE

Enables controlled mobilization of elbow.
medial and lateral epicondyle fractures.
After Surgeries.
elbow dislocation,elbow arthrolysis.



SU-3007

2 Size (S-M)

2 Side (Right-Left)

ELBOW IMMOBILIZER



Made of sponge that is covered with the terry fabric from the inside.
Stabilize the Elbow thanks to the metal pieces in it

USAGE

Protections at impacts.
Fractures.
Protection during training and sport activities.

SU-3008

Size:

UNISIZE

WRIST AND FINGER SPLINTS



WRIST SPLINT LONG

Made of elastic cotton fabric specially designed to constrain the movements of the wrist thanks to the aluminum rod inside.

Provides better stabilization thanks to the elastic bandage on the wrist.

USAGE

Provides immobilization and constrains to the movements of the wrist.

Arthritis.

Rheumatism.



SU-4001

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25 cm



WRIST SPLINT (NEOPRENE) LONG



Made of neoprene fabric that provides warmth and protection.

Constrains movements of the wrist thanks to the aluminum support inside.

USAGE

Wrist tunnel cases.

Tendonitis, tenosynovitis.

After taking off gypsum.

After Surgeries.

SU-4002

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25 cm



WRIST SPLINT (NEOPRENE) Right and Left

Made of neoprene fabric that provides warmth and protection.

Constrains movements of the wrist thanks to the aluminum support inside.

USAGE

Wrist tunnel cases.

Tendonitis, tenosynovitis.

After taking off gypsum.

After Surgeries.



SU-4004

Size:

UNISIZE

WRIST SPLINT WITH THUMB SUPPORT LONG



Made of elastic cotton fabric specially designed to constrain the movements of the wrist thanks to the aluminum rod inside.

The thumb supporting rod stabilizes it in the rest position.

Provides better stabilization thanks to the elastic bandage on the wrist.

USAGE

Protection and stabilization of the wrist and the thumb.

Arthritis and rheumatism.

Wrist tunnel cases.

Pains in the thumb and wrist

SU-4005

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25 cm



WRIST SPLINT (NEOPRENE) With Thumb Support Long

Made of neoprene fabric that provides warmth and protection. It has an aluminum support with a bandage to spread the pressure all over the hand.

USAGE

Stabilizing the thumb.
Wrist tunnel cases.
Tendonitis, tenosynovitis.
After taking off gypsum.
After Surgeries.



SU-4006

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25cm



WRIST SPLINT SHORT



Made of elastic cotton fabric specially designed to constrains the movements of the wrist thanks to the aluminum rod inside.

USAGE

Protection and stabilization of the wrist .
Arthritis and rheumatism.
Wrist tunnel cases.
Pains in wrist

SU-4007

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25cm



WRIST SPLINT SHORT NEOPRENE

Made of neoprene fabric that provides warmth and protection. It has an aluminum support to spread the pressure all over the hand.

USAGE

Protection and stabilization of the wrist .
Arthritis and rheumatism.
Wrist tunnel cases.
Pains in wrist



SU-4008

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25cm



WRIST SPLINT WITH THUMB SUPPORT (Short) Neoprene Fabric



Stabilizing the thumb's joint.
Thumb ligamen's problems.

USAGE

Stabilizing the thumb's joint.
Thumb ligamen's problems.

SU-4009

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25cm



ADJUSTABLE WRIST SPLINT (NEOPRENE)

It is manufactured from perforated fabric. It is a type of splint that provides immobilization of the wrist in neutral position. It has an aluminium support inside and 2 plastic supports at the top. There is a bandage that provides stabilization in the laced system around the wrist and can be used in both hands.

USAGE

Carpal Tunnel Syndrome
In tendon and soft tissue injuries
Rehabilitation after surgical intervention



SU-4011

Size:
UNISIZE

NEOPRENE WRIST BANDAGE



It is a simple wrist bandage. It provides protection and heat to the wrist in case of wrist injuries, sprains and trauma. It improves blood circulation by maintaining the body temperature and thus speeds up the healing process.

USAGE

In wrist strains
To prevent sports injuries
In carpal dislocation

SU-4013

Size:
UNISIZE

NEOPRENE HAND WRIST SUPPORT

It is a simple wrist bandage that put thumb to inside. It provides protection and heat to the wrist in case of wrist injuries, sprains and trauma. It improves blood circulation by maintaining the body temperature and thus speeds up the healing process.

USAGE

Twists, Trumas, bruises.
Protection while doing sports.
Protection at work when the jobs involves using the wrist a lot.



SU-4014

Size:
UNISIZE

THUMB SUPPORT (DE QUEUVAIN)



It is a type of splint that immobilize thumb. The plastic atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

USAGE

CMC and MCP joint irritation
In de quervain syndrome
In ligament damage

SU-4015

Size:
UNISIZE

THUMB SUPPORT

Made of thermoplastic material, It is a type of splint that immobilize thumb. The plastic atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

USAGE

CMC and MCP joint irritation
In de quervain syndrome
In ligament damage



SU-4016

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm



WRIST BAND WITH THUMB SUPPORT



It is a type of splint that immobilize thumb. The aluminium atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

USAGE

CMC and MCP joint irritation
In de quervain syndrome
In ligament damage

SU-4017

Size:

UNISIZE

WRIST BAND WITH THUMB SUPPORT LONG

It is a type of splint that immobilize thumb. The aluminium atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

USAGE

CMC and MCP joint irritation
In de quervain syndrome
In ligament damage



SU-4018

Size:

UNISIZE

WRIST SPLINT WITH DOUBLE SIDE



There is an anatomically shaped aluminium support that provides immobilization of the wrist in the neutral position. There is a bandage that increases the stabilization around the wrist. It is a standard size and can be personally adjusted with velcro.

USAGE

Carpal Tunnel Syndrome
In tendon and soft tissue injuries
Rehabilitation after surgical intervention

SU-4019

Size:

UNISIZE

WRIST SPLINT WITH THUMB SUPPORT NEOPRENE

There is an aluminium support under the wrist and near the thumb. It can use to both hands by changing aluminium support under the wrist. There is a velcro that increases the stabilization around the wrist.

USAGE

In de quervain syndrome
Carpal Tunnel Syndrome
Rehabilitation after surgical intervention
In wrist and thumb calcification
In tendon and soft tissue injuries



SU-4020

Size:
UNISIZE

THERMOPLASTIC RELIEVING SPLINT



Hand and wrist orthosis manufactured from thermoplastic material. Keeps the fingers and the wrist in recovery position. Enables immobilization of wrist, used for resting position following plaster treatment of fractures and painful positions of the wrist.

USAGE

For bone fractures of the wrist after plaster casting
For resting in soft tissue injuries
Carpal Tunnel Syndrome
Rehabilitation after surgical intervention
In tendon and soft tissue injuries

SU-4021

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm



THERMOPLASTIC ANTI-SPASTICITY HAND SPLINT

Made of thermoplastic covered by terry fabric from the inside. Stabilizes the wrist and the fingers in a relaxing opened position.

USAGE

Spasticity
Wrist problems.



SU-4022

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm



THERMOPLASTIC ANTI-SPASTICITY HAND SPLINT LUX



Made of thermoplastic covered by terry fabric from the inside. Stabilizes the wrist and the fingers in a relaxing opened position.

USAGE

Spasticity
Wrist problems.

SU-4023

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm



THERMOPLASTIC DYNAMIC HAND-FINGER SPLINT

Made of thermoplastic covered by terry fabric from the inside.
Pulls up the fingers and let them move a bit to help training them.

USAGE
Rehabilitation.
Extensor tendon.s.



SU-4024

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm



THERMOPLASTIC KLEINERT HAND-FINGER SPLINT



Made of thermoplastic covered by terry fabric from the inside.
Pulls up the fingers and let them move a bit to help training them.

USAGE
Rehabilitation.
Flexure tendon.

SU-4025

Size:


S	15-17 cm
M	17-19 cm



MALLET FINGER SPLINT

Made of high quality plastic.
Stabilize the upper part of the finger and prevent it from moving.

USAGE
Mallet finger deformities.
Interruption of tendons.
Fractures and bruises.



SU-5001

Size:

S	4.5 - 5 cm
M	5 - 6 cm
L	6 - 7 cm
XL	7 - 8 cm
2XL	8 - 9 cm

FROG SPLINT



Made of the Soft aluminum witch has the ability to reshape as wished.
The inner side is covered with the plastoz witch helps stabilizing the finger.

USAGE
Finger fractures.
Strings and nerve injuries.
Finger injuries and pain.

SU-5002

Size:

S	4 - 6 cm
M	6 - 8 cm
L	8 - 10 cm

BASEBALL SPLINT

Made of the Soft aluminium witch has the ability to reshape as wished.
The inner side is covered with the plastoz witch helps stabilizing the finger.

USAGE

Finger fractures.
Strings and nerve injuries.
Finger injuries and pain.



SU-5003

Size:

S 4 - 6 cm

M 6 - 8 cm

L 8 - 10 cm

BUTONIER FINGER SPLINT

Made of the Soft aluminium witch has the ability to reshap as wished.
The inner side is covered with the plastoz witch helps stabilizing the finger.

USAGE

Rehabilitation.
Strings injuries.
Finger movement training.



SU-5004

Size:

S 4 - 6 cm

M 6 - 8 cm

L 8 - 10 cm

STRAIGHT FINGER SPLINT

Made of the Soft aluminium which has the ability to reshape as wished.
The inner side is covered with the plastoz witch helps stabilizing the finger.

USAGE

Finger fractures.
Strings and nerve injuries.
Finger injuries and pain.



SU-5005

Size:

Standard

CORSETS



ELASTIC LUMBOSTAD CORSE-WITH BELT (26 CM)

Made of elastic cotton fabric specially designed to reduce sweating.
 Has 2 metal and 2 plastic supporters reduces the pressure on the discs.
 Has elastic belt provides extra stabilization.
 It's 26 cm High.
 The metal supporters can be removed while washing.
USAGE
 Lumbar problems.
 Chronic back pain.



SU-6001

Size:

- S 70-80 cm
- M 80-90 cm
- L 90-100 cm
- XL 100-110 cm
- XXL 110-120cm



LUMBOSACRAL CORSET 26 CM




Made of elastic cotton fabric specially designed to reduce sweating.
 Has 2 metal and 2 plastic supporters reduces the pressure on the discs.
 It's 26 cm High.
 The metal supporters can be removed while washing.
USAGE
 Sliding discs.
 Chronic back pain.
 Befor and after surgeries.
 Lumbar spine pain.

SU-6002

Size:

- S 70-80 cm
- M 80-90 cm
- L 90-100 cm
- XL 100-110 cm
- XXL 110-120cm



ELASTIC LUMBOSTAD CORSE-WITH BELT (32 CM)


Made of elastic cotton fabric specially designed to reduce sweating. Has 2 metals and 2 plastics supporters reduces the pressure on the discs. Has elastic belt provides extra stabilization. It's 32 cm High.
 The metal supporters can be removed while washing.
USAGE
 Sliding lumbar and sacral discs.
 Chronic back pain. Lumbar problems.
 Befor and after surgeries.
 Reduces the pressure on the discs.



SU-6003

Size:

- S 70-80 cm
- M 80-90 cm
- L 90-100 cm
- XL 100-110 cm
- XXL 110-120cm



DORSO LUMBAR CORSET (THORACO LUMBAR)




Made of elastic cotton fabric specially designed to reduce sweating.
 Has 2 metal and 2 plastic supporters at the lower back.
 Supports the spine and the shoulders thanks to its design.
 Comfortable during daily use.
USAGE
 Spine problems.
 Osteoporosis.
 The small humping.
 Support before and after Surgeries.

SU-6004

Size:

- S 70-80 cm
- M 80-90 cm
- L 90-100 cm
- XL 100-110 cm
- XXL 110-120cm



CORSETS

POSTUREX BANDAGE

Made of elastic cotton fabric specially designed to reduce the Sweating.
Has a plastic supporting reduces the ability to bend.
Has side laces helps pulling the shoulders to keep the body's stabilization.
It can be easily worn under the clothes.
USAGE
Shoulder and back muscles pain.
Other muscle problems.
Long computer and car driving hours.
Back bags carrying



SU-6005

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



MAGNETIC BACK CORSET




Adjustable straps design to fit all sizes the flexibility that gives you the freedom of movement. It has a magnetic belt that supports the spine to provide stability for those suffering from lower back pain

USAGE
Great for lower back pain and supporting the back. It also keeps lower back in an upright or extended position it reduces pain in muscles, ligaments, discs and nerve roots of the lower back

SU-6006

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



ABDOMINAL CORSET


Made of elastic cotton fabric specially designed to reduce sweating.
Stabilize the abdomen.
Can be adjusted as needed.
Reduces the the pains in the abdomen area.
Designed to fit under clothes.
USAGE
Post abdominal surgery and birth
Peritoneal injury
Abdominal sagging
Post umbilical and ventral hernia surgery



SU-6007

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



CHEST CORSET



Made of elastic cotton fabric specially designed to reduce sweating.
Its inside is covered by smooth non-allergenic fabric.
Designed to fit under clothes.

USAGE
Sternal Support after coronary surgery
Costa and sternum fractures
Post chest surgery and trauma

SU-6008

Size:

S	80-90 cm
M	90-100 cm
L	100-110 cm
XL	110-120 cm
XXL	120-130cm



NEOPRENE LUMBOSACRAL CORSET WITH BELT

Made of neoprene fabric that provides warmth and protection.
 Has 2 metals and 2 plastics supporters reduces the pressure on the discs.
 Has a belt provide extra warmth and pressure.
USAGE
 Sliding lumbar and sacral discs.
 Back pain.
 Prevent edemas.
 Enhance the blood circulation.



SU-6009

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



NEOPRENE LUMBOSACRAL CORSET 26 CM




Made of neoprene fabric that provides warmth and protection. Has 2 metals and 2 plastics supporters reduces the pressure on the discs. It's 26 cm High.
 Has a belt provide extra warmth and pressure.
USAGE
 Sliding discs.
 Chronic back pain.
 Befor and after surgeries.
 Lumbar spine pain.

SU-6010

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



ABDOMINAL CORSET (NEOPRENE)

Made of neoprene fabric that provides warmth and protection.


USAGE
 Provides warmth to the abdomen and the lower part of the back.
 Helps burning the fats and losing weight during doing sports.



SU-6011

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



HYPEREXTENSION CORSET




Made of a special designed aluminum. The inner side is covered with a smooth non-allergic sponge fabric. The 3 points system stabilize the spine and prevent its movement which speeds up the healing process.

USAGE
 Stable compression fracture of lumbar and lower thoracic vertebrae
 After inter-vertebral disc surgery
 After non-stable vertebral body fracture surgery. Non-operable spinal tumor and its metastases. Post-operation orthopedic care and orthosis
 After spinal decompression surgery

SU-6012

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



HYPEREXTENSION CORSET ADJUSTABLE

Made of a special designed aluminum. The inner side is covered with a smooth non-allergic sponge fabric. The special design helps modifying it as needed. Stabilizes the spine and keeps it straight.

USAGE

Stable compression fracture of lumbar and lower thoracic vertebrae
After inter-vertebral disc surgery
After non-stable vertebral body fracture surgery. Non-operable spinal tumor and its metastases. Post-operation orthopedic care and orthosis
After spinal decompression surgery



SU-6013

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



COLOSTOMY CORSET



Made of elastic cotton fabric specially designed to reduce sweating. Has a side hole in it so the colostomy base can be attached. Designed to fit under clothes.


USAGE

Protecting the ostomy collection bag
Port-ostomy rehabilitation and support purposes
Decrease risks of parastomal hernia

SU-6014

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



ABDOMINAL CORSET (WOOL)

Made of 2 layer wool fabric.

USAGE


Provides warmth to the abdomen and the lower part of the back. cal and ventral hernia surgery



SU-6015

Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
XXL	110-120cm



WOOL CORSET (OPEN FRONT SIDE)



Made of 2 layer wool fabric..


USAGE

Provides warmth to the abdomen and the lower part of the back.

SU-6016

Size:

S	80-90 cm
M	90-100 cm
L	100-110 cm
XL	110-120 cm
XXL	120-130cm



INGUINAL RAPTURE BELT RIGHT OR LEFT

Made of a special cotton fabric.
Does not effect the daily life.
Provides support after surgeries.
Fits perfectly under clothes.

USAGE
Single sided inguinal hernia



SU-6018

Size:
Standard

INGUINAL RAPTURE BELT (DOUBLE)

Made of a special cotton fabric.
Does not effect the daily life.
Provides support after surgeries.
Fits perfectly under clothes.

USAGE
Double sided inguinal hernia



SU-6019

Size:
Standard

(NEOPRENE) SLIMMING SHORT

Made of neoprene fabric
Contains pores that help for ventilation.
effectively improve blood circulation, Keeps the heat, reduces the pressure, reduces the severity of muscle's pain. Helps to slim the hip, thigh and torso. Designed for comfort, support and protection for hamstring and thigh injury. Zipper provide ease at wear.

USAGE
Providing heat and compression
Inguinal problems
Slimming around the waist and hip line
Injuries to the thigh



SU-6020

Size:
S 70-80 cm
M 80-90 cm
L 90-100 cm
XL 100-110 cm
XXL 110-120cm



KNEE SUPPORTS



(NEOPRENE) KNEE SUPPORT CLOSED PATELLA

Made of neoprene fabric
Contains pores that help for ventilation.
Support knee joint, reduce pain by relieving the pressure on the joint of the knee.
Sustains body's own release of heat effectively improve blood circulation, promote metabolism, enhance immunity also applies compression, prevents edema.

USAGE

Arthritis. Intra articular liquid increase
Pain induced by knee cap
Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities
When heat and compression is desired



SU-7001

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



(NEOPRENE) KNEE SUPPORT CLOSED WITH SILICONE PAD



Made of neoprene fabric
Contains pores that help for ventilation.
Support knee joint, reduce pain by relieving the pressure on the joint of the knee.
The perimeter of the patella is closed and supported with a silicone pad.
Sustains body's own release of heat effectively improve blood circulation, promote metabolism, enhance immunity also applies compression, prevents edema.

USAGE

Arthritis. Intra articular liquid increase
Pain induced by knee cap
Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities
When heat and compression is desired

SU-7002

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



(NEOPRENE) KNEE SUPPORT PATELLA SUPPORT

Made of neoprene fabric
Contains pores that help for ventilation.
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Can be used daily as it does not compress your knee cap.
Recommended for people who have pain in their knee cap with and difficulty when climbing stairs or squatting.

USAGE

Arthritis. Intra articular liquid increase
Pain induced by knee cap
Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities
When heat and compression is desired



SU-7003

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



(NEOPRENE) KNEE SUPPORT PATELLA WITH LIGAMENT SUPPORT- SHORT



Made of neoprene fabric
Contains pores that help for ventilation.
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Flexible wire on the sides increase stabilization and decrease the pressure on the knee.

USAGE

Arthritis, bursitis, tendinitis
Sprains and stains
Meniscus problems
Post-op or trauma for protection purposes
Ligament instability
Sporting injuries

SU-7004

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



(NEOPRENE) KNEE SUPPORT PATELLA WITH LIGAMENT SUPPORT- LONG

Made of neoprene fabric
Contains pores that help for ventilation.
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Flexible wire on the sides increase stabilization and decrease the pressure on the knee.

USAGE

Arthritis, bursitis, tendinitis
Sprains and stains. Meniscus problems
Post-op or trauma for protection purposes
Ligament instability
Sporting injuries



SU-7005

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



(NEOPRENE) KNEE SUPPORT HALF OPENED PATELLA With LIGAMENT SUPPORT



Made of neoprene fabric
Contains pores that help for ventilation.
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Flexible wire on the sides increase stabilization and decrease the pressure on the knee.

USAGE

Arthritis, bursitis, tendinitis
Sprains and stains. Meniscus problems
Post-op or trauma for protection purposes
Ligament instability. Sporting injuries

SU-7006

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



KNEE SUPPORT WITH HINGED STABILIZING JOINTS

Made of neoprene fabric
Contains pores that help for ventilation.
Aluminum joints on both sides provide maximum support, allow controlling the motion and prevent knee hyperextension. Velcro fasteners on the upper and lower calves increase stabilization.

USAGE

Ligament injury and instability
Meniscus lesions
1 and 2 degree sprains and stains
Post-op or trauma for protection purposes



SU-7007

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



HALF OPENED KNEE SUPPORT WITH HINGED STABILIZING JOINTS



Made of neoprene fabric
Contains pores that help for ventilation.
Aluminum joints on both sides provide maximum support, allow controlling the motion and prevent knee hyperextension. Velcro fasteners on the upper and lower calves increase stabilization.

USAGE

Ligament injury and instability
Meniscus lesions
1 and 2 degree sprains and stains
Post-op or trauma for protection purposes

SU-7008

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



KNEE SUPPORT WITH HINGED STABILIZING JOINTS- FULL OPEN FRONT

Made of neoprene fabric
Contains pores that help for ventilation.
Aluminum joints on both sides provide maximum support, allow controlling the motion and prevent knee hyperextension.
Velcro fasteners on the upper and lower calves increase stabilization.
Velcro fasteners on the upper and lower calves increase stabilization.
Since it can be completely open from the front side, provides ease of use to the elderly, heavy weight and post-op patients.

USAGE

Ligament injury and instability
Meniscus lesions
1 ve 2 degree sprains and stains
Post-op or trauma for protection purposes



SU-7009

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



KNEE SUPPORT PATELLA WITH CROSSED LIGAMENT SUPPORT-NEOPRENE



Made of neoprene fabric
Contains pores that help for ventilation.
Contains opened the patella supported by a pad. Medial (inside) and lateral (outside) Flexible wire offer overall knee joint stability.

Transverse straps adhering with 3-point principle provide stabilization for ant. Cruciate ligaments and patella.

USAGE

Ant. Cruciate ligament injury and knee instability
Patellar stabilization
Meniscus lesions
Post-op or trauma for protection and rehabilitation purposes

SU-7010

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



KNEE SUPPORT FOR PATELLAR TENDON

Made of neoprene fabric that provides warmth and protection.
The special pad supports the patellar tendon and with the applied compression reduces pain.

USAGE

Os-good Schlatter's disease
Patellar tendinitis
Patellar chondromalacia
Patellofemoral pain



SU-7011

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



PATELLAR TENDON BAND UNISIZE



Made of neoprene fabric that provides warmth and protection.
The special pad supports the patellar tendon and with the applied compression reduces pain.

Does not restrict joint motion.
Providing comfortable all-day wear.
The knee support strap is made with neoprene so it provides a comfortable and flexible fit.
Used for both knees.

USAGE

Os-good Schlatter's disease
Patellar tendinitis
Patellar chondromalacia
Patellofemoral pain

SU-7012

**Size:
Standard**

WOOL KNEE SUPPORT

Manufactured of wool for warmth and compression

USAGE

Knee pain
Liming
Rheumatic pains
Temperature for



SU-7013

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



KNEE SUPPORT OPEN PATELLA



Made of neoprene fabric
Contains pores that help for ventilation.
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Can be used daily as it does not compress your knee cap.
Recommended for people who have pain in their knee cap with and difficulty when climbing stairs or squatting.

USAGE

Arthritis. Intra articular liquid increase
Pain induced by knee cap
Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities. When heat and compression is desiredg injuries

SU-7014

Size:
Standard

KNEE SUPPORT OPEN PATELLA WITH LIGAMENT

Made of neoprene fabric
Contains pores that help for ventilation.
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Can be used daily as it does not compress your knee cap.
Recommended for people who have pain in their knee cap with and difficulty when climbing stairs or squatting.

USAGE

Arthritis. Intra articular liquid increase
Pain induced by knee cap
Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities. When heat and compression is desiredg injuries



SU-7015

Size:
Standard

ADJUSTABLE ANGLE KNEE ORTHESIS WITH BAR LUX



Use for protected activity during post-op rehab following ligamentous surgical repairs and reconstruction. Controlled R.O.M. for post-traumatic, post-surgical repair of human knee. It can be extended and shortened with the buttons on it.

USAGE

Controlled and gradual joint mobilization after ligament injury operations
Meniscal repairs. For knee rest and stabilize at the desired angle after surgery.
In knee recurvatum injuries
Stable femoral fractures

SU-7020

Size:
Standard

ADJUSTABLE ANGLE KNEE ORTHESIS WITH BAR

Adjustable-degree knee orthosis, manufactured from aluminium bar and sponge laminated cloth with a terry inner side. Flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) and extension (-15°-0°, 15°, 30°, 45°, 60°, 75°, 90°) angle adjustment possible. Provides exceptional support and stability following knee injury. Use for protected activity during post-op rehab following ligamentous surgical repairs and reconstruction. Controlled R.O.M. for post-traumatic, post-surgical repair of human knee.

USAGE

Controlled and gradual joint mobilization after ligament injury operations. Meniscal repairs
For knee rest and stabilize at the desired angle after surgery. In knee recurvatum injuries
Stable femoral fractures



SU-7021

Size:
Standard (47/58cm)

FUNCTIONAL KNEE BRACE

It is a functional rigid brace that is made of light-weight aluminium material designed according to the anatomical knee shape, providing flexion and extension limitation. Flexion (0°, 15°, 30°, 60°, 90°, 120°) and extension (0°, 15°, 30°) degrees can be adjusted. It provides the most favorable stabilization with the principle of 4 points. It prevents forward and back tibial displacement and prevents knee hyperextension..

USAGE

Instability caused by collateral ligament rupture
Conservative treatment of combined instabilities
instabilities
Isolated rupture of ant. cruciate ligament or combined ant. and post. Cruciate ligament ruptures



SU-7022

Size:
Standard

NEOPRENE KNEE SUPPORT WITH CONTROLLED MECHANISM

Neoprene fabric which Provides firm support and promotes flexibility and it also offers therapeutic heat to soothe. aching joints. Controls extension and flexion degrees of the knee movement. Screws that are adjusted in Flexion in degree. (0-15-30-60-90-120)and Extension (0-5-15-30) provides external. support by protecting knee junction. Better control, comfort, and increased compliance.

And reduces the pain by relieving the load on the joint.

USAGE

Ligament injuries and instability,. Meniscus lesions, 1st and 2nd degree sprains and strains.
Post-op or trauma for protection and rehabilitation purposes. Degenerative joint disease where extension/flexion. control is required.



SU-7023

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



KNEE IMMOBILIZER

Manufactured of sponge laminated with a terry cloth inner coating, This orthosis immobilizes the leg at 20°. Aluminum supports provide immobilization in extension, stabilization in the posterior, medial and lateral. Knee immobilizers are ideal for knee injuries and as a postoperative aid as they are easily removed. Standard size.

Lengths available at 45cm-55cm and 65 cm.

Thanks to the three-piece pad construction, this orthosis can be easily adjusted to any leg circumference.

USAGE

Recovery after ligament injury

Meniscus Post-op recovery

Knee joint in trauma cases

Osteochondral fractures

Acute sprains/strains of the knee

Patellar dislocations

Post-op recovery, protection and immobilization



SU-7024

Size:
45 CM



SU-7025

Size:
55 CM



SU-7026

Size:
65 CM

THIGH SUPPORT

Made of neoprene fabric
Contains pores that help for ventilation.
It supports thigh's muscles
effectively improve blood circulation,
Keeps the heat, reduces the pressure,
reduces the severity of muscle's pain

USAGE

Muscle injuries to the thigh
Thigh pain
Protection during and after the sportive activities
When heat and compression is needed



SU-7027

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



CALF SUPPORT



Made of néoprène fabric
Contains pores that help for ventilation.
It supports thigh's muscles
effectively improve blood circulation,
Keeps the heat, reduces the pressure, re-
duces the severity of muscle's pain

USAGE

Muscle injuries to the calf
Calf pain
Protection during and after the sportive activities
When heat and compression is needed

SU-7028

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



HIP SUPPORT

It is manufactured to support thigh. It re-
duces pain by radiate heat to thigh area.
Velcros on it provides easy usage. Include
belt and so it provides comfortable usage
to patients for daily usage.

USAGE

Muscle injuries to the thigh
Thigh pain
Protection during and after the sportive ac-
tivities
When heat and compression is needed



SU-7029

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



HIP ABDUCTION ORTHESIS



Hip abduction Orthosis anatomically de-
signed and manufactured of thermoplastic
material, providing hip stabilization.
Parts in contact with the skin are coated
with terry cloth. Flexion adjustment possi-
ble at every 15 degrees, starting from 15
up to 120 degrees.

Abduction angle is adjustable.

USAGE

Stabilization of total hip prosthesis or revi-
sion
Femoral head resection
Femoral or acetabular osteotomy
Recurrent hip dislocation

SU-7030

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



ANKLE SUPPORTS



ANKLE SUPPORT (CROSS STRAP)

Made of neoprene fabric
Contains pores that help for ventilation.
The figure 8 elastic transverse strap provides extra stabilization and compression.
The soft neoprene blend retains body heat to improve circulation and promote healing process.

USAGE

1st degree twist and tonus
Minor fractures
Compression on edema
Protection during and after the sportive activities



SU-8001

Size:

S	24-26 cm
M	26-28 cm
L	28-31 cm
XL	31-34 cm
XXL	34-37cm



ANKLE SUPPORT WITH SILICONE PAD (MALLEOLUS)



Made of neoprene fabric
Contains pores that help for ventilation.
The Crossed lacing provides extra pressure and stabilization thanks to the rubber pads on the sides.

USAGE

Ligament instability and rapture
Edema prevention
Post-op or trauma for protection and rehabilitation purposes

SU-8002

Size:

S	24-26 cm
M	26-28 cm
L	28-31 cm
XL	31-34 cm
XXL	34-37cm



LIGAMENT ANKEL SUPPORT

Made of neoprene fabric
Contains pores that help for ventilation.
The aluminium around the ankle provides high leveled stabilization.
The aluminium can be modified as wanted.

The Crossed lacing provides extra pressure and stabilization.

USAGE

Ligament instability and rapture
Edema prevention
Post-op or trauma for protection and rehabilitation purposes



SU-8003

Size:

S	24-26 cm
M	26-28 cm
L	28-31 cm
XL	31-34 cm
XXL	34-37cm



FIGURE OF EIGHT BANDAGE



Made of neoprene fabric
Contains pores that help for ventilation.
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Flexible wire on the sides increase stabilization and decrease the pressure on the knee.

USAGE

Arthritis, bursitis, tendinitis
Sprains and stains
Meniscus problems
Post-op or trauma for protection purposes
Ligament instability
Sporting injuries

SU-8004

Size:
Standard

ANKLE SUPPORT (PLASTIC + LACED)

Made of neoprene fabric
Contains pores that help for ventilation.
Thanks to two specially designed plastic plates located on the malleolus, it provides strong support to the ankle.

USAGE

Provides stronger stabilization to the ankle with serious ankle sprains.
Also effective with habitual ankle sprains, chronic ankle pains and can be used as a protector with contact sports.



SU-8005

Size:

S	24-26 cm
M	26-28 cm
L	28-31 cm
XL	31-34 cm
XXL	34-37cm



ANKLE SUPPORT (CROSS STRAP)



Made of neoprene fabric
Contains pores that help for ventilation.
The figure 8 elastic transverse strap provides extra stabilization and compression.
The soft neoprene blend retains body heat to improve circulation and promote healing process.

USAGE

1st degree twist and tonus
Minor fractures
Compression on edema
Protection during and after the sportive activities

SU-8006

Size:
Standard

ANKLE BRACE WITH AIR PAD

Made with specially designed plastic molds.
Thanks to the air pads inside, it fits to the foot anatomically and provides maximum stabilization.

USAGE

Functions as a stabilizer in cases of serious ankle sprains and chronic ankle pains and as a protector in contact sports.
Also applicable in cases of ankle fractures and in the period after plastering process.



SU-8007

Size:
Standard

ANKLE BRACE WITH PLASTIC PAD



The hard outer shell has been made with specially designed plastic molds.
Thanks to the soft pads inside, it fits to the foot anatomically and provides maximum stabilization.

USAGE

With its rigid exterior, provides good stabilization in pronation and supination actions.
Provides full support for ankle ligament injuries.

SU-8008

Size:
Standard

ROM WALKER

It comes with 3 pads (12.5°) so it can be adjusted as needed.

Physiological walk due to special sole.

Plantar flexion: 0°, 10°, 20°, 30°, 40°

Dorsiflexion: 20°, 10°, 0°

USAGE

After ligament, soft tissue and tendon surgeries.

For post operative and follow up treatments of the talus, calcaneus and medial malleolus fractures.

For early functional therapy after conservative or surgical treatment for rupture for the Achilles tendon



SU-8009

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



ACHILL BOOT



It comes with 3 pads (12.5°) so it can be adjusted as needed.

Physiological walk due to special sole.

USAGE

After ligament, soft tissue and tendon surgeries.

For post operative and follow up treatments of the talus, calcaneus and medial malleolus fractures.

For early functional therapy after conservative or surgical treatment for rupture for the Achilles tendon

SU-8010

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



REFLEX AFO

Ankle orthosis made of plastic material.

Provides static dorsiflexion assistance and lateral stability for entire ankle/foot area.

Provides ease of use with velcro fasteners going around the ankle.

Use a scissor to trim insole for exact shoe fit.

May be modified with a heat gun.

USAGE

Foot-drop deformity caused by nerve injuries such as hemiplegia, lumbar disc hernia

When 90 degrees ankle angle is needed



SU-8011

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



FOOT GUARD



T shaped ankle orthosis made of thermoplastic material with terry cloth coating.

Ensures the fix 90 degrees stand, prevent rotation of the ankle with T Shape.

USAGE

Stabilizing the position of the leg on hospital patients

SU-8012

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



HALLUKS VALGUS NIGHT SPLINT

Foam laminated fabrics, the inner surface is made of terry cloth. Position of the thumb to the desired position is hidden inside the fabric with aluminum splint can be adjusted. designed to be used at Night (in sleep) .

USAGE
Helps to keep the big toes aligned, designed to use during night time.



SU-8013

Size:
Standard

HALLUKS VALGUS DAY SPLINT

Velcro and elastic bandage is made of neo-prene fabric. Day (shoes, slippers ... in) designed to be used at night.

USAGE
Helps to keep the big toes aligned, designed to use during day time.



SU-8014

Size:
Standard

CASTHOE

with adjustable Velcro tapes can also be applied to both feet.

USAGE
Used instead of shoes in case of plaster feet



SU-8015

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



FOOT CARE



SILICON INSOLE

made of high quality of silicone material. Specially designed to fit figure of the foot. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf.

Distributes the standing body weight evenly and reduces load and stress on the joints. Absorbs the reaction force of the floor.

USAGE

Pes planus
Plantar fasciitis
Heel pain
Metatarsal pain



SU-8017

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



SILICONE INSOLE WITH ARCH SUPPORT

Specially designed to fit figure of the foot. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Arc support supplement for low base and related pain relief.

Distributes the standing body weight evenly and reduces load and stress on the joints.

USAGE

Provides anatomical support to the sole and enables even distribution of pressure on the sole.
Has a relaxing and comforting effect.



SU-8021

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



PERFORATED SILICONE INSOLE WITH ARCH SUPPORT

Specially designed to fit figure of the foot. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Arc support supplement for low base and related pain relief. Distributes the standing body weight evenly and reduces load and stress on the joints.

USAGE

Provides anatomical support to the sole and enables even distribution of pressure on the sole. Has a relaxing and comforting effect.



SU-8022

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



COVERED SILICONE INSOLES



Specially designed to fit figure of the foot and covered with fabric. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Distributes the standing body weight evenly and reduces load and stress on the joints. Absorbs the reaction force of the floor.

USAGE

Pes planus. Plantar fasciitis
Heel pain. Metatarsal pain

SU-8023

Size:

S	34-37 cm
M	38-39 cm
L	40-41 cm



CORK INSOLES

Longitudinal and metatarsal arch support manufactured from mixture of special cork and rubber material, covered with leather surface. Distributes the standing body weight evenly and reduces load and stress on the joints.



SU-8032

Size:

From 35
To 44



ARC REINFORCED RUBBER INSOLE



Specially designed to fit figure of the foot and made of rubber material. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Arc support supplement for low base and related pain relief. Distributes the standing body weight evenly and reduces load and stress on the joints.

SU-8033

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



CORK INSOLE WITH ARCH SUPPORT 3/4 LENGTH

Longitudinal and metatarsal arch support manufactured from mixture of special cork and rubber material. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Arc support supplement for low base and related pain relief. Distributes the standing body weight evenly and reduces load and stress on the joints.



SU-8034

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



STEEL INSOLE



made of stout leather and supported with stainless steel bridge has a longitudinal and metatarsal support.

Distributes the standing body weight evenly and reduces load and stress on the joints.

Rigid surface support does not allow the longitudinal arc to collapse.

USAGE

Pes planus-planovalgus
Pes cavus
Traverse arc support

SU-8035

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



SILICONE METATARSAL PAD WITH RING

It is made of fabric material that does not harm the skin. The position of the silicone ring is adjusted by inserting it on the 2nd finger. Metatarsals relieve pain by removing the pressure on their heads.

USAGE

Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.



SU-8118

Size :

Standard

3/4 INSOLES (SILICONE)



Middle degree of hardness, is made from 100% silicon. 3/4 in size and offers comfortable use.

USAGE

Provides anatomical support to the sole and enables even distribution of pressure on the sole. Has a relaxing and comforting effect.

SU-8018

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm



SILICON HEEL CUP

made of high quality of silicone material. Specially designed to fit figure of the foot. Softened calcaneal epin area lessens the load on the heels which causes pain. Lessens the pain at foot and calf.

USAGE

Heel spur
Plantar fasciitis
Avescular necrosis of calcaneus (sever's disease)
Heel pain



SU-8019

Size:

S	34-37 cm
M	38-39 cm
L	40-41 cm



SILICON TOE SEPERATOR



made of high quality of silicone material. Comfortable for daily use. Prevents deviation of the toe to the lateral in hallux valgus deformity.

USAGE

Helps proper positioning of toes in the shoes and prevents overlapping.

SU-8020

Size :

Standard

SILICON WEDGE

It is made of silicone material and the upper surface is covered with fabric. In accordance with the problem, it is placed on the inner or outer side of the shoe. Due to its adhesive nature, it can be fixed within the shoe.

USAGE

- 1-Calcaneus varus,
- 2-Calcaneus valgus,
- 3-O-Leg, X-Leg,
- 4-Meniscus degeneration.vity.



SU-8024

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



SILICONE HEEL CUP FABRIC COATED WITH EPINE



Manufactured from silicon material, covered with fabric. Softened calcaneal epin area prevents the load pressing onto a single point. Due to its adhesive nature, it can be fixed within the shoe.

USAGE

- 1-Heel spur,
- 2-Plantar fasciitis,
- 3-Avascular necrosis of calcaneus (sever's disease),
- 4-Heel pain.

SU-8025

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



SILICONE METATARSAL PAD

It is made of silicone material and the upper surface is covered with fabric. It is placed within the proximal of the metatarsus heads. It relieves the pressure of the metatarsus heads and reduces pain. Due to its adhesive nature, it can be fixed within the shoe.

USAGE

- 1-Over-sensitivity of metatarsus head,
- 2-Metatarsalgia,
- 3-Callus developing under the metatarsus,
- 4-Morton's foot- hallux rigidus-claw toe deformity,
- 5-Aceptycal necrosis of metatarsus head (Freiberg's disease).



SU-8026

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm



SILICON ARC SUPPORT



It is made of silicone material and the upper surface is covered with fabric. It supports the longitudinal arch and reduces the increased tension in the sole. It starts from the calcaneus' inner tubercle and continues until the rear of the metatarsus heads. Due to its adhesive nature, it can be fixed within the shoe.

USAGE

- 1-Pes planus,
- 2-Pes cavus,
- 3-When the longitudinal arch needs to be supported.

SU-8027

Size:

S	34-37 cm
M	38-39 cm
L	40-41 cm



SILICONE HAMMER TOE CUSHION (PAIR)

It is made of silicone material that does not harm the skin. It is used by passing on hammer fingers. It helps to eliminate deformation by applying support over the joint.

USAGE

- Over-sensitivity of metatarsus head
Metatarsalgia
Callus developing under the metatarsus



SU-8028

Standard Size :

Left / Right

SLICONE METATARSAL PAD WITH RING

Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.

USAGE

Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.



SU-8016

Size :
Standard

SLICONE METATARSAL PAD WITH TOESEPERATOR



Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.

USAGE

Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.

SU-8029

Size :
Standard

SILICON BUNION SHIELD

Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.

USAGE

Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.



SU-8030

Size :
Standard

SILICONE TOE CAP




made of high quality of silicone material. Comfortable for daily use. Prevents deviation of the toe to the lateral in hallux valgus deformity.

USAGE

Is made of silicone. Best to use in the day-time. Used to heal deviation in lateral toe Hallux Valgus.

SU-8031

Size:



VARICOSE SOCKS CALF (opened or closed toe)

USAGE

Varicosis Stockings apply a decreasing compression towards the calf and upper leg whose external pressure is highest on the ankle.

Thus, they expedite the blood flow by providing a regular pressure distribution in legs and feet.

Patellar dislocations

Post-op recovery, protection and immobilization

CLASS 1

Knitted elegantly and helps to relax
Pressure offers of 14 mmHg

SU-8037

CLASS 2

Knitted elegantly and helps to relax
Pressure offers of (20-30 mmHg)


SU-8040

CLASS 3

Knitted elegantly and helps to relax
Pressure offers of (30-40 mmHg)

SU-8043

Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
cB	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
cB1	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
cC	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
cD	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
cE	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
cF	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
cG	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92



VARICOSE SOCKS THIGH (opened or closed toe)

USAGE

Varicosis Stockings apply a decreasing compression towards the calf and upper leg whose external pressure is highest on the ankle.

Thus, they expedite the blood flow by providing a regular pressure distribution in legs and feet.

Patellar dislocations

Post-op recovery, protection and immobilization

CLASS 1

Knitted elegantly and helps to relax
Pressure offers of 14 mmHg

SU-8038

CLASS 2

Knitted elegantly and helps to relax
Pressure offers of (20-30 mmHg)


SU-8041

CLASS 3

Knitted elegantly and helps to relax
Pressure offers of (30-40 mmHg)

SU-8044

Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
cB	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
cB1	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
cC	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
cD	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
cE	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
cF	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
cG	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92



VARICOSE SOCKS (PANTY)

USAGE

Varicosis Stockings apply a decreasing compression towards the calf and upper leg whose external pressure is highest on the ankle.

Thus, they expedite the blood flow by providing a regular pressure distribution in legs and feet.

Patellar dislocations

Post-op recovery, protection and immobilization

CLASS 1

Knitted elegantly and helps to relax
Pressure offers of 14 mmHg

SU-8039

CLASS 2

Knitted elegantly and helps to relax
Pressure offers of (20-30 mmHg)


SU-8042

CLASS 3

Knitted elegantly and helps to relax
Pressure offers of (30-40 mmHg)

SU-8045

Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
cB	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
cB1	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
cC	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
cD	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
cE	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
cF	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
cG	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92




ANTIEMBOLISM SOCKS

USAGE

prevent leg swelling and helps with leg discomfort, improves blood circulation, commonly known as blood clots, advised for post surgery or prolonged bedrest.

ANTIEMBOLISM SOCKS CALF

Prevent varicose veins
It accelerates blood circulation. It increases the blood flow of the veins



SU-8046

ANTIEMBOLISM SOCKS THIGH

Prevent varicose veins
It accelerates blood circulation. It increases the blood flow of the veins



SU-8047


ANTIEMBOLISM SOCKS THIGH WITH BELT

Muscles pain.
Neck arthritis.
Sliding discs.
After the traumas and bruises



SU-8048

Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
cB	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
cB1	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
cC	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
cD	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
cE	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
cF	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
cG	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92





**PEDIATRIC
PRODUCTS**

SOFT CERVICAL COLLAR

It provides superior comfort and support to the neck through anatomical designed. It refrains from excessive cervical activity and provides extra warmth around the neck. Compatibility MR.

USAGE

Muscles pain.
Neck arthritis.
Sliding discs.
After the traumas and bruises



SU-9001

Size :
Standard

CERVICAL COLLAR (NELSON) FOR KIDS



It provides superior comfort and support to the neck through anatomical designed. It helps the patient during the treatment by restricting the movements in the neck area. Compatibility MR.

USAGE

Muscles pain.
Neck arthritis.
Sliding discs.
After the traumas and bruises

SU-9002

Size :
Standard

CERVICAL COLLAR (PHILEDELPHIA) FOR KIDS

It is a 2-piece neck collar. It helps the patient during the treatment by limiting flexion, lateral and rotational movements of the neck. It has a chin support. Compatibility MR.

USAGE

Soft tissue trauma.
Ligaments and bones injuries.
gliding discs.
Stabilizations in the cases of the unconsciousness.



SU-9003

Size :
Standard

CERVICAL COLLAR (PHILEDELPHIA - TRACHEOSTOMY) FOR KIDS



It is a 2-piece neck collar. It helps the patient during the treatment by limiting flexion, lateral and rotational movements of the neck. It is a collar with trachea open. Compatibility MR.

USAGE

Soft tissue trauma.
Ligaments and bones injuries.
Sliding discs.
Stabilizations in the cases of the unconsciousness.

SU-9004

Size :
Standard

CERVICAL COLLAR WITH CHIN SUPPORT (PLASTIC) FOR KIDS

Made of PVC with chin support in a design that can be modified as wanted to stabilize the neck and reduce the pressure on the spine

USAGE

Muscles pain.
Neck arthritis.
Sliding discs.
Cases of radiculopathy
After the traumas and bruises



SU-9005

Size:
Standard

CERVICAL COLLAR (PLASTIC) FOR KIDS



Made of PVC in a design that can be modified as wanted to stabilize the neck and reduce the pressure on the spine

USAGE

Sliding discs
Neck arthritis
After the traumas and bruises

SU-9006

Size:
Standard

ARM SLING MESH FABRIC FOR KIDS

manufactured from perforated fabric. Product's design provides patient's arm to be fixed parallel to the near trunk. It can be used for both sides.

USAGE

Fractures and bruises of the arm.
Bruised shoulder.



SU-9007

Size:

S	16-22 cm
M	22-28 cm
L	28-34 cm



CLAVICULA BANDAGE KIDS



Made of comfortable cloth covered from the inside with a layer of soft cloth which does not cause allergies, helps to stabilize the wishbone thanks to belts that pull back the shoulders to the extent required to give the appropriate positioning.


USAGE

Clavicle dislocations
Clavicle fractures

SU-9008

Size:

S	40-50 cm
M	50-60 cm
L	60-70 cm



ELBOW SUPPORT KIDS

Made of neoprene fabric that provides warmth and protection. Thanks to the three-headed silicone pad inside, it applies constant pressure and compression to the epycondil heads.

USAGE

Lateral epicondylitis.
Internal articular fluid.
Elbow pain



SU-9009

Size:

S	16-22 cm
M	22-28 cm
L	28-34 cm



WRIST BANDAGE KIDS



manufactured of neoprene fabric
Provides protection for the wrist, sprains and bruises wrist
Retains body heat to improve bloody circulation and accelerate wound healing

USAGE

Twists, Trumas, bruises.
Protection while doing sports.
Protection at work when the jobs involves using the wrist a lot.

SU-9010

Size :

Standard

WRIST SPLINT KIDS

Made of neoprene fabric that provides warmth and protection. Constrains movements of the wrist thanks to the aluminium support inside.

USAGE

Wrist tunnel cases.
Tendonitis, tenosynovitis.
After taking off gypsum.
After Surgeries.



SU-9011

Size:

S	9-11cm
M	11-13 cm
L	13-15 cm



WRIST SPLINT WITH THUMB SUPPORT (KIDS)



There is an aluminium support under the wrist and near the thumb. It can use to both hands by changing aliminium support under the wrist. There is a velcro that increases the stabilization around the wrist.

USAGE

Wrist tunnel cases.
Tendonitis, tenosynovitis.
After taking off gypsum.
After Surgeries.

SU-9012

Size:

S	9-11cm
M	11-13 cm
L	13-15 cm



NEOPREN KNEE SUPPORT WITH PAD FOR KIDS

Made of neoprene fabric that provides warmth and protection. Used during doing sports.

USAGE

Supporting the Knee and protecting it from injuries.



SU-9013

Size:

S	25-28 cm
M	28-31 cm
L	31-34 cm



NEOPREN KNEE SUPPORT WITH LIGAMENTS (KIDS)



Made of neoprene fabric that provides warmth and protection. The flexible underwire on the sides increase stabilization and thus decrease the load on the knee.

The soft neoprene blend retains body heat to improve circulation and promote healing.

USAGE

Arthritis, bursitis. Increase of intra articular liquid. Pain induced by knee cap
Meniscus problems. 1ST degree twist and ligament instability. Post-op or trauma for

SU-9014

Size:

S	25-28 cm
M	28-31 cm
L	31-34 cm



PEDIATRIC KNEE IMMOBILIZER (KIDS)

It provides immobilization in extension. It is designed anatomically and includes aluminum supports to provides stabilization in the posterior, medial and lateral. Thanks to the three-piece pad design, it can be easily adjusted for leg circumference.

USAGE

Recovery after ligament injury. Meniscus Post-op recovery. Knee joint in trauma cases
Ostrochondral fractures. Acute prains/strains of the knee. Patellar dislocations. Post-op recovery, protection and immobilization



SU-9015

Size:

S	25-28 cm
M	28-31 cm
L	31-34 cm



ANKLE SUPPORT KIDS



Made of neoprene fabric that provides warmth and protection.. Velcro fasteners around the ankle provide ease at wear. Supports the ankle thus reduces the pain. The soft neoprene blend retains body heat to improve circulation and promote healing process.

USAGE

1st degree twist and tonus
Minor fractures
Compression on edema
Protection during and after the sportive activities

SU-9016

Size:

S	18-20 cm
M	20-22 cm
L	22-24 cm



FIGURE OF EIGHT BANDAGE

Made of neoprene fabric
Contains pores that help for ventilation.
The Crossed lacing provides extra pressure and stabilization thanks to the rubber pads on the sides.

USAGE

Ligament instability and rapture
Edema prevention
Post-op or trauma for protection and rehabilitation purposes



SU-9017

Size :
Standard

PAVLIK SUPPORT



Made of 3 cm cotton belts, non-elastic, stabilize the hip and the foot in a perfect position.

USAGE

Congenital Dislocation

SU-9018

Size:

S 0-3 Month

M 3-6 Month

L 6-9 Month

FREJKA ORTHOSIS

Manufactured from 12 mm top quality molded plastazote foam.
Holes allow air to circulate, soothing support and relief.
Latex-free, non-toxic and hypoallergenic plastazote foam material will not cause skin irritation. Velcro straps are used for easy donning and doffing. Plastazote foam material also is radiolucent which allows X-ray, it possible to use CT, MRI. Positions at the hip joint in flexion and abduction.

USAGE

Hip subluxations



SU-9019

Size:

S 0-3 Month

M 3-6 Month

L 6-9 Month

UMBILICAL HERNIA TRUSS



made of flexible cotton fabric. The pad applies pressure to the inguinal hernia. Perineal strap prevents sliding. It comes in standard size but adjustable.

USAGE

UMBILICAL HERNIA

SU-9020

Size :
Standard

INGUINAL RAPTURE BELT

Made of a special cotton fabric.
Does not effect the daily life.
Provides support after surgeries.
Fits perfectly under clothes.

USAGE
Single sided inguinal hernia



SU-9021

Standard Size :
Left/Right

INGUINAL RAPTURE PEDIATRIC BELT (DOUBLE)



Made of a special cotton fabric.
Does not effect the daily life.
Provides support after surgeries.
Fits perfectly under clothes.

USAGE
Double sided inguinal hernia

SU-9022

Size :
Standard

DENNIS BROWN (MOVABLE/STABLE)

The child can do independent movement of the legs such as kicking or crawling, decrease the risk of the boot sliding off of the foot and regional rash. The child will be active and sleeps more comfortably.

USAGE
Used for hip dislocation and pes equinovarus deformity.



SU-9023

Size :
Standard

DENNIS BROWN SHOES



SU-9024

Size:
From 17
to 27



**OTHER
PRODUCTS**



SILICON BREAST

It is used to provide a natural appearance when chest tissue is removed after mastectomy, to reduce pain caused by muscle imbalances, to protect the surgical wound from external bumps and crushes. The laundry is made of cotton and non-allergenic special fabric and silicone material. Shoulder straps are softened for comfortable use.


The pockets on both sides allow the prosthesis to be used as right and left.



SU-100

Size:
Standard

SILICON BREAST'S BRA



The prosthetic bra, which is very comfortable and stylish to wear, is specially designed for breast prosthesis. It supports the prosthesis and weight and provides a sense of confidence and comfortable use. It is made of cotton fabric, the inner surface of which does not irritate the skin. Shoulder straps are softened for comfortable use. The pockets on both sides allow the prosthesis to be used both left and right.

SU-101

Size:
Standard

PORTABLE CLOSED SEAT LIFTER (WITHOUT COVER)

Made of 3 cm cotton belts, non-elastic, stabilize the hip and the foot in a perfect position.

USAGE
Congenital Dislocation



SU-102

Size:
Standard

PORTABLE CLOSED SEAT LIFTER (WITH COVER)



Latex-free, non-toxic and hypoallergenic plastizode foam material will not cause skin irritation. Velcro straps are used for easy donning and doffing. Plastizode foam material also is radiolucent which allows X-ray, it possible to use CT, MRI. Positions at the hip joint in flexion and abduction.

USAGE
Hip subluxations

SU-103

Size:
Standard

SILICON EXERCISE BALLS

Made of high quality silicone material which does not cause allergies and has no smell. Designed to withstand the pressure applied on, without losing their shape.

Available in three different colors

USAGE

- Used so as to rehabilitation.
- Three different hardnesses existing.
- Red – Hard, Blue – Medium, Green – Soft



SU-104

Size:
Standard

VISCO NECK PILLOW



Size :
60*40*11/8 cm

SU-200

VISCO NECK PILLOW



Size :
60*37*11/9 cm

SU-201

VISCO CLASSIC NECK PILLOW



Size :
60*40*12 cm

SU-202

VISCO CLASSIC NECK PILLOW



Size :
60*40*15 cm

SU-203

SNORING REDUCING VISCO PILLOW



Size :
55*40*12/10 cm

SU-204

APNEA PILLOW



Size :
55*40*12/10 cm

SU-205

VISCO NECK PILLOW WITH GEL



Size :
57*37*11/7 cm

SU-206

VISCO CLASSIC NECK PILLOW WITH GEL



Size :
60*40*15 cm

SU-207

PILLOWS

VISCO NECK PILLOW WITH GEL



Size :
60*37*11/9 cm

SU-208

FLAT CLASSIC PILLOW WITH GEL



Size :
55*40*11 cm

SU-209

VISCO TRAVEL PILLOW



Size :
31*27*7 cm

SU-210

VISCO TRAVEL PILLOW



Size :
31*27*10 cm

SU-211

VISCO SEAT CUSHION



Size :
43*9 cm

SU-212


VISCO SEAT CUSHION (OPEN)



Size :
43*43*9,5 cm

SU-213


VISCO SQUARE CUSHION



Size :
45*43*7 cm

SU-214

VISCO U SQUARE CUSHION



Size :
41*39*6 cm

SU-215

VISCO ANATOMIC SEAT



Size :
42*40*9.5 cm

SU-216

VISCO COCCYX CUSHION



Size :
44*36*9.5 cm

SU-217


VISCO POSTURE SEATING CUSHION



Size :
44*38*13.7 cm

SU-218


VISCO SQUARE CUSHION WITH GEL



Size :
45*43*7 cm

SU-219

VISCO ANATOMIC SEAT WITH GEL



Size :
42*40*9.5 cm

SU-220

BETWEEN LEG VISCO PILLOW



Size :
25*19*14 cm

SU-221

BACK AND WAIST SUPPORT PILLOW



Size :
54*42*9 cm

SU-222

VISCO LUMBAR CUSHION



Size :
38*32*12 cm

SU-223

VISCO LUMBAR CUSHION



Size :
35*34*12 cm

SU-224


VISCO LUMBAR CUSHION



Size :
34*32*8 cm

SU-225

VISCO LUMBAR CUSHION



Size :
39*34*9 cm

SU-226

VISCO REFLUX PILLOW



Size :
65*60*24 cm

SU-227

Partner[®]
Let Us Improve Your Health

**İKİTELLİ OSB MAH. 6. CAD.
BEYAZ TOWER No:1 Office
No:118 BAŞAKŞEHİR
İSTANBUL - TURKEY**

**+90 553 588 1688
+90 212 671 1917
www.partner-medical.com
Info@Partner-Medical.Com**