

# Partner<sup>®</sup>

*Let Us Improve Your Health*

**PRODUCTS CATALOGUE  
2019 - 2020**



**Chairy<sup>®</sup>**

**Walky<sup>®</sup>**

**Bubbly<sup>®</sup>**

**Body<sup>®</sup>**

# Company Profile

PARTNER, Mobility Aids and orthopedics designer & manufacturer since 2014.

Our Vision is to be a major player in the world of medical equipment. The PARTNER has specialized over the years in the manufacturing and designing of orthopedics, wheelchairs, Commode chairs, comfort chairs, Crutches, Canes, Rollator and any other product that may accompany a patient or user in maintaining its autonomy and mobility.

Our European development and design teams are equipped with precision high-tech tools resulting in reliability and quality for each product we produce. Today, PARTNER is a group of more than 100 employees with subsidiaries in Europe and Asia and delivers products sold in more than 15 countries worldwide.



**TURKEY**

# Certifications

**CERTIFICATE OF REGISTRATION**



SYNDICATE OF INTERNATIONAL SYSTEM

This Certificate has been awarded to

**PARTNER TRADING DIŞ TİCARET  
LİMİTED ŞİRKETİ**

İKİTELLİ OSB MAH. 6. CAD. BEYAZ TOWER NO:1/118  
BAŞAKŞEHİR / İSTANBUL / TÜRKİYE

In recognition of the organization's Managements System which complies with

**ISO 9001:2015**

The scope of activities covered by this certificate is defined below  
DESIGN, SUBCONTRACT MANUFACTURER, SALES, STORAGE, DISTRIBUTION,  
IMPORT & EXPORT OF MOBILITY AIDS PRODUCTS FOR DISABLED PEOPLE AND  
ORTHOPEDIC PROSTHESES, ORTHOSIS, SUPPORTS AND BRACES PRODUCTS,  
ACCESSORIES AND PRODUCTION OF ORTHOPEDIC, REHABILITATION  
DEVICES AND MOBILITY AIDS PRODUCTS

Certificate Number: SISTURQ1020191456  
Date of Issue of Original Certificate: 14.10.2019  
Date of Issue of Latest Certificate: 14.10.2019  
Expiry Date: 13.10.2020



Managing Director

Note : This certificate is valid only if produced with the continuation letter after the surveillance is carried out successfully.

The Organization's documentation and Implementation has been reviewed and found to comply with the relevant standard rules. This certificate of Registration is based on the evaluation of the mentioned scope given above. Organization is responsible for maintaining the responsibilities of the relevant standard rules. Any significant changes in the scope of the certification or standard referred above render this certificate invalid. This is an accredited certificate issued by SIS Certifications Pvt. Ltd. sanctioned for issue by International Accreditation Services, 3060 Saturn Street Suite 100 Brea, California 92821-1732, USA.

Email us : support@siscertifications.com, info@siscertifications.co.in, Call: +90-9654721648  
Web: - http://www.siscertifications.co.in, www.siscertifications.com  
The status of this certificate can be verified on "http://www.siscertifications.co.in".

Issue No.: 01






## Certificate

### CE

#### CE DECLARATION OF CONFORMITY

Manufacture Name: PARTNER TRADING DIŞ TİCARET LİMİTED ŞİRKETİ

Address: İKİTELLİ OSB MAH. 6. CAD. BEYAZ TOWER NO:1/118 BAŞAKŞEHİR / İSTANBUL / TÜRKİYE

The undersigned company certifies under its sole responsibility that the item groups specified below satisfies the requirements of the MEDICAL DEVICES 93/42/EEC which apply to it. The item groups identified below has been subject to internal manufacturing checks with monitoring of the final assessment by PARTNER TRADING DIŞ TİCARET LİMİTED ŞİRKETİ

**Description of Product:**  
ORTHOPEDIC PRODUCTS

**Product Model:**  
CERVICAL COLLARS, KNEE BRACES, ARM SLING, CORSETS, ELBOWS, HAND, WRIST SPLINTS, HIP LEG BRACES, ANKLE BRACES, FOOT SUPPORTS, SILICON PRODUCTS, VARICOSE STOCKINGS, ELASTIC BANDAGES, KNITTED PRODUCTS, ACCESSORIES AND OTHER ORTHOPEDICS PRODUCTS

**Related Directives:**  
MEDICAL DEVICES 93/42/EEC

**Regulations Applied acc. To Harmonized Standards:**  
EN ISO 15223-1:2013 SYMBOLS TO BE USED WITH MEDICAL DEVICE LABELS, LABELLING AND INFORMATION TO BE SUPPLIED - GENERAL REQUIREMENTS  
EN 1041:2008+A1:2013 INFORMATION SUPPLIED BY THE MANUFACTURER OF MEDICAL DEVICES  
EN ISO 14971:2013 MEDICAL DEVICES - APPLICATION OF RISK MANAGEMENT TO MEDICAL DEVICES

Certificate Number : CE/19596  
Certificate Issue Date : 14.10.2019  
Certificate Validity : 13.10.2020





QCS International Certifications Services Pvt. Ltd  
Certificate of validity information E-mail: info@qcsint.com  
Head to contact



## Certificate

### CE

#### CE DECLARATION OF CONFORMITY

Manufacture Name: PARTNER TRADING DIŞ TİCARET LİMİTED ŞİRKETİ

Address: İKİTELLİ OSB MAH. 6. CAD. BEYAZ TOWER NO:1/118 BAŞAKŞEHİR / İSTANBUL / TÜRKİYE

The undersigned company certifies under its sole responsibility that the item groups specified below satisfies the requirements of the MEDICAL DEVICES 93/42/EEC which apply to it. The item groups identified below has been subject to internal manufacturing checks with monitoring of the final assessment by PARTNER TRADING DIŞ TİCARET LİMİTED ŞİRKETİ

**Description of Product:**  
MOBILITY AIDS, ASSISTIVE DEVICES AND REHABILITATION DEVICE

**Product Model:**  
MANUAL WHEELCHAIRS, ELECTRIC WHEELCHAIRS, COMMODOE CHAIRS, CANES, CRUTCHES, WALKERS, HOSPITAL BEDS, HOSPITAL FURNITURE, ANTI BED SORES MATTRESS, PHYSIOTHERAPY AND REHABILITATION DEVICES AND OTHER MOBILITY AIDS PRODUCTS, DEVICES AND ACCESSORIES

**Related Directives:**  
MEDICAL DEVICES 93/42/EEC

**Regulations Applied acc. To Harmonized Standards:**  
EN ISO 15223-1:2013 SYMBOLS TO BE USED WITH MEDICAL DEVICE LABELS, LABELLING AND INFORMATION TO BE SUPPLIED - GENERAL REQUIREMENTS  
EN 1041:2008+A1:2013 INFORMATION SUPPLIED BY THE MANUFACTURER OF MEDICAL DEVICES  
EN ISO 14971:2013 MEDICAL DEVICES - APPLICATION OF RISK MANAGEMENT TO MEDICAL DEVICES

Certificate Number : CE/19625  
Certificate Issue Date : 14.10.2019  
Certificate Validity : 13.10.2020





QCS International Certifications Services Pvt. Ltd  
Certificate of validity information E-mail: info@qcsint.com  
Head to contact

## CERTIFICATE of Registration



This is to Certify that the  
Medical Devices -- Quality Management System  
of

**PARTNER TRADING DIŞ TİCARET LİMİTED ŞİRKETİ**

İKİTELLİ OSB MAH. 6. CAD. BEYAZ TOWER NO:1/118  
BAŞAKŞEHİR / İSTANBUL / TÜRKİYE

has been independently assessed and is compliant  
with the requirements of

### ISO 13485:2016

This Certificate is applicable to the following product or service ranges:  
DESIGN, SUBCONTRACT MANUFACTURER, SALES, STORAGE, DISTRIBUTION,  
IMPORT & EXPORT OF MOBILITY AIDS PRODUCTS FOR DISABLED PEOPLE AND  
ORTHOPEDIC PROSTHESES, ORTHOSIS, SUPPORTS AND BRACES PRODUCTS,  
ACCESSORIES AND PRODUCTION OF ORTHOPEDIC, REHABILITATION  
DEVICES AND MOBILITY AIDS PRODUCTS

Certificate No.: TR51295H

Date of initial registration	12 October 2019
Date of this Certificate	12 October 2019
Surveillance audit on or before	11 October 2020
Recertification Due / Certificate expiry	11 October 2022

This Certificate is property of Staunchly Management & System Services Pvt. Ltd. and remains valid subject to satisfactory surveillance audits.



Director

STAUNCHLY MANAGEMENT & SYSTEM SERVICES LTD.  
Suite 48, 88-90 Hutton Garden, London, EC1N 8PN.

Phone: +44 20 800 0199  
Email: info@staunchlyservices.com Web: www.staunchlyservices.com  
For precise and updated information concerning the present certificate mail to info@staunchlyservices.com  
This Certificate is the property of Staunchly Management & System Services Private Limited and shall be returned immediately when demanded

SMS/190A/17/REV02

# Index Contents

## Wheelchairs

- Manual Wheelchairs ..... 5-12
- Electrical Wheelchairs ..... 13-15
- Commode Wheelchairs ..... 16-18

**Bath Bench** ..... 20

**Walking Aids** ..... 22-25

**Bubble mattress** ..... 27

**Electric Hospital Bed** ..... 28



# Partner<sup>®</sup>

*Let Us Improve Your Health*

[www.partner-medical.com](http://www.partner-medical.com)



**WheelChair**  

---

.....

## P - 100 A








# Chairy®

### Product features

- Steel powder coated construction.
- Sit back cushion red Oxford cloth.
- Fixed armrest and backrest tube .
- Fixed footrest with plastic footplate.
- 8" (20cm) front castors.
- 22" (56cm) rear spoke wheels with solid tires.
- Double brake.
- Foldable wheelchair.



Note : All dimensin in ( cm , kg )

Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 45cm - 18"	 40	 88	 63	 92	 15.2 / 13.5	 115

## P - 100

# Chairy®








### Product features

- Frame Steel, single cross brace.
- Surface treatment powder coating.
- Front wheel 8"x1" front castor.
- 24" 9 spoke MAG wheel with grey PU tire.
- Padded flip back armrest with desk length.
- Swing away detachable footrest.
- Front fork Steel.
- Height adjustable Foot plate Plastic.
- Nylon seat and backrest.
- Two leg straps for safety.
- Foldable wheelchair.



**Optional:**

- Seat Belt - Elevating Legrests - Pneumatic Rear Wheels - Seat Cushion - Different Colors / Sizes

Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 40cm - 16" 46cm - 18" 50cm - 20"	 40	 92	 66	 112	 19.8/17.8	 115

Note : All dimensin in ( cm , kg )

## P - 101

# Chairy®








### Product features

- Powder coating grey wheelchair frame with single cross braces.
- 24" (61cm) 9 spoke MAG wheel with grey PU tire.
- 8" x 2" (20cm x 5cm) front castor with grey PU tires.
- Padded flip back armrest with desk length.
- Swing away wide pace footrest (nib spring style).
- Leg support. • Plastic hand rim.
- Powder coating steel front fork.
- Plastic footplate.
- Nylon seat and backrest.
- Foldable wheelchair.



#### Optional:

- Seat Belt - Elevating Legrests - Pneumatic Rear Wheels - Seat Cushion - Different Colors / Sizes

Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 40cm - 16" 46cm - 18" 50cm - 20"	 40	 91	 66	 112	 19.8/17.8	 115

Note : All dimensin in ( cm , kg )



## P - 100 C

# Chairy®

### Product features

- Steel pediatric wheelchair.
- Yellow powder coating wheelchair frame.
- 20" spoke wheel with pneumatic tires.
- 6" front castor with PU tires.
- Powder coating steel hand rim.
- Padded flip back up armrest with desk length.
- Swing away detachable blue footrest with leg strip.
- Black powder coating aluminum alloy front fork.
- Height adjustable Plastic footplate.
- Nylon seat and backrest with sponge.
- Against backward overturn device.
- Foldable wheelchair.



**Seat width**

**Seat depth**

**Wheelchair height**

**Wheelchair width**

**Wheelchair length**

**Weight**

**Maximum capacity**

GW / NW



35cm - 14"

35

94

59

91

15/20

75

Note : All dimensin in ( cm , kg )

## P - 102

## Chairy®








### Product features

- Steel wheelchair.
- Powder coating wheelchair frame with single cross brace.
- Height adjustable 24" (61cm) 9 spoke MAG wheel with grey PU tire.
- 8x1" (20cm x 2.5cm) front castor with black solid rubber tires.
- Black plastic hand rim.
- Detachable full-length armrest.
- Vinyl upholstery with pocket & flame-resistant seat.
- Swing away Nylon footrest.
- Height adjustable alloy front fork.
- Leg support.
- Foldable wheelchair.



#### Optional:

- Seat Belt - Elevating Legrests - Pneumatic Rear Wheels - Seat Cushion - Different Colors / Sizes

Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 46cm - 18" 50cm - 20"	 40	 90	 66	 112	 21/18.5	 115

Note : All dimensin in ( cm , kg )

## P - 103








### Product features

- Powder coated grey Aluminum frame with single cross braces.
- 24" (61cm) 9-spoke MAG wheel with grey PU tire.
- 8x2" (20cm x 5cm) grey front castor with PU tires.
- Plastic front fork.
- Black plastic hand rim.
- Slope armrest with soft PU arm pad.
- Nylon seat with sponge.
- Swing in and away footrest.
- Height adjustable plastic footplate.
- Foldable wheelchair.



#### Optional:

- Seat Belt
- Pneumatic Rear Wheels
- Seat Cushion
- Different Colors / Sizes

Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 40cm - 16" 46cm - 18" 50cm - 20"	 46	 91	 61	 103	 21.5/19.0	 120

Note : All dimensin in ( cm , kg )

## P - 104 R

# Chairy®








### Product features

- Steel Wheelchair.
- Powder coating wheelchair frame with double braces.
- 24" spoke rear wheels with grey PU tires.
- 8"x2" front castors with grey PU tires.
- Powder coating steel hand rim.
- Padded detachable armrest with full length.
- Elevating footrest with padded largest.
- Powder coating Steel front fork.
- Height adjustable plastic footplate.
- Vinyl seat and backrest.
- Reclining backrest.
- Foldable wheelchair.



#### Optional:

- Seat Belt
- Pneumatic Rear Wheels
- Seat Cushion
- Different Colors / Sizes

Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 46cm - 18" 50cm - 20"	 50	 119	 65	 130	 29.7/27.7	 120

Note : All dimensin in ( cm , kg )

## P - 105CP

# Chairy<sup>®</sup>

### Product features

- Ultra-light steel frame.
- Powder coating wheelchair frame.
- Easy to adjust chest holders to fix the user body upright position.
- Reclining high back and adjustable headrest.
- Stylish mag style wheels light up this product.
- The frame allows tilting between 3°-22° up to user demand.
- 16" spoke rear wheels.
- 6"x2" front castors.
- Powder coating steel front fork.
- Detachable armrest.
- Elevating footrest.
- Foldable wheelchair.



**Seat width**



46 cm - 18"

**Seat depth**



40

**Wheelchair height**



103

**Wheelchair width**



61

**Wheelchair length**



103

**Weight**

GW / NW



22.2 / 20.4

**Maximum capacity**



120

Note : All dimensin in ( cm , kg )

## P - 340w

# Chairy®

Enjoy the independence of a power wheelchair with the portability of a lightweight manual wheelchair.

### Product features

Frame	Steel, powder coating
Controller	MAIKONG 50A
Motor	340 x 2 LINIX
Battery	Led - acid 2x12 36 Ah
Charger	HP 24V, 4A
Max. speed	7 Km/h
Travel milage	28 Km
Gradeability	10 degree
Armrest	Flip back slope armrest
Seat	Detachable cushion
Footrest	Swing away , heel strap
Other	Anti-tip, backrest half folding
Front wheel	8" PU caster
Rear wheel	12.5" Penumatic tire



**Seat width**



43cm - 17"  
46cm - 18"  
48cm - 19"

**Seat depth**



43

**Wheelchair height**



97

**Wheelchair width**



63

**Wheelchair length**



108

**Weight**

GW / NW



60/56

**Maximum capacity**



136

Note : All dimensin in ( cm , kg )

## P - 250w








# Chairy®

Enjoy the ultimate and durable driving experience with the P-250w most advanced with featured MAIKONG Control technology and innovative controls, easy to fold wheelchair.

### Product features

Frame	Steel, powder coating
Controller	WUYANG 45A
Motor	250W x 2
Battery	2x12, 20 Ah
Charger	AC220 50~60Hz
Max. speed	6 Km/h
Travel milage	15 Km
Gradeability	6 degree
Armrest	Flip back slope armrest
Footrest	Swing in and away
Seat width	17" / 18" / 19"
Front wheel	10"
Rear wheel	18"
Side panel	Plastic
Other	With anti-tip



Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 43cm - 17" 46cm - 18" 48cm - 19"	 43	 97	 45	 78	 60/56	 120

Note : All dimensin in ( cm , kg )

## P - 250wa

# Chairy<sup>®</sup>

Made from durable, lightweight steel making it extremely durable and able to hold up to 130 kg. The simple folding design allows for convenient storage and easy transport.

### Product features

Frame	Steel, powder coating
Controller	MAIKONG
Motor	250W x 2
Battery	2x12, 20 Ah
Charger	AC220 50~60Hz
Max. speed	6 Km/h
Travel milage	20 Km
Footrest	Swing in and away
Gradeability	5 degree
Armrest	Flip back slope armrest
Front wheel	8" PU folaming solid
Rear wheel	16" Rubber inflatable
Side panel	Plastic



Seat width



43cm - 17"

Seat depth



43

Wheelchair height



95

Wheelchair width



65

Wheelchair length



100

Weight  
GW / NW



40.8/35.7

Maximum capacity



130

Note : All dimensin in ( cm , kg )



## CO - 1

### Product features

- Steel powder painting frame.
- United attendant brakes for convenient attendant usage.
- commode padded seat
- Front wheel 8"x1" front PVC solid castor.
- 23" Solid PVC rear wheels.
- Vinyl seat and backrest.
- Leg strap for safety.
- Foldable wheelchair.



**Seat width**

**Seat depth**

**Wheelchair height**

**Wheelchair width**

**Wheelchair length**

**Weight**

**Maximum capacity**



45 cm - 18"



43



90



66



102



7.75 / 15.6



130

Note : All dimensin in ( cm , kg )

**CO - 2**

**CO - 3**



**Product features**

- Foldable Chromed steel frame.
- Plastic commode seat.
- Fixed armrest and elevating footrest.
- Solid castors and rear wheels.
- Reclining high back with handle under the armrest.

**Product features**

- Foldable Chromed Steel frame.
- Plastic commode seat.
- Flip-up armrest and detachable footrest.
- Solid castors and solid rear wheels.

Note : All dimensin in ( cm , kg )

	Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
CO-2	44 cm - 13"	43	129	64	117	28.7 / 25.2	100
CO-3	45 cm - 18"	43	90	60	106	21.4 / 19.6	100

CO - 4

CO - 5

















Product features

- Steel frame ( Chrome paltded ).
- leather upholstery.
- Flip down armrest.
- Detachable footrest.
- Paking brake mechanism.
- Plastic commode seat.

Product features

- Aluminum Frame.
- PE seat and back.
- round bucket.
- 3" wheels with parking mechanism.

Note : All dimensin in ( cm , kg )

	Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight	Maximum capacity
CO-4	 44 cm - 17"	 43	 92	 55	 82	 8.8	 136
CO-5	 44 cm - 17"	 45	 89	 55	 84	 5	 136

**Chairy<sup>®</sup>**

**B - 1**



**B - 2**



**B - 3**



**B - 4**



	Seat width	Seat depth	Chair height	Chair width	Chair length	Weight GW / NW	Maximum capacity
<b>B-1</b>	50 - 20"	30	78 - 84	40	40	12.8 / 10.8	100
<b>B-2</b>	40 cm - 16"	33	72 - 85	40	44	8.2 / 7.6	100
<b>B-3</b>	50 cm - 20"	30	39 - 49	40	40	4.4 / 3.6	100
<b>B-4</b>	84 - 130   33" - 51"	36	80 - 90	45	38	10.2 / 9.3	100

Note : All dimensin in ( cm , kg )

The logo features the word "Walky" in a bold, dark blue, sans-serif font. A thin red curved line is positioned below the text. The entire logo is centered within a white circle that has a thin teal border. A vertical white line with a circular dot at the bottom is located on the left side of the page.

**Walky<sup>®</sup>**

**Walking Aids**  

---

## WA - 1

Walky®

**Product features**

- Made of Aluminum.
- Reflective elements ensure safe movement in dark area.
- Height adjustable.

**Use:**

Walking aid which improves stability and maintains upright ambulation,  
especially useful when an injury or disability

Height (cm)	79-94
Weight (kg)	0.3
Maximum capacity (kg)	115



## WA - 2

**Product features**

- Made of Aluminum.
- Reflective elements ensure safe movement in dark area.
- Height adjustable.

**Use:**

Walking aid which improves stability and maintains upright ambulation,  
especially useful when an injury or disability

Height (cm)	75-98
Weight (kg)	0.56
Maximum capacity (kg)	115



## WA - 3

Walky®

**Product features**

- Made of Aluminum.
- Reflective elements ensure safe movement in dark area.
- Height adjustable.

**Use:**

Walking aid which improves stability and maintains upright ambulation, especially useful when an injury or disability.

Height (cm)	67.5-90
Weight (kg)	0.36
Maximum capacity (kg)	115



## WA - 4

**Product features**

- Made of Steel.
- Reflective elements ensure safe movement in dark area.
- Height adjustable.

**Use:**

Walking aid which improves stability and maintains upright ambulation, especially useful when an injury or disability.

Height (cm)	72-95
Weight (kg)	0.77
Maximum capacity (kg)	120



WA - 5

Walky®

Product features

- Made of Aluminum.
- Height adjustable. • Packed in pairs
- Armpit and hand parts covered in polyurethane foam, warm in contact.

Use:

Walking aid which improves stability and maintains upright ambulation.

Size	S	M	L
Height (cm)	84-114	114-134	134-154
Weight (kg)	1.6	1.7	1.8
Maximum capacity (kg)	115	115	115



WA - 6

Product features

- 3 Functions walker can be used as reciprocal, fixed or wheeled walker.
- Foldable light Aluminum frame.
- Stable and durable construction.
- Removable rubber walker tips.
- Additional removable front tips with 12 cm diameter.
- Ergonomic handles. • Height adjustable.

Use:

Walking aid which improves stability and helps maintains upright ambulation.

Height (cm)	77-96
Weight (kg)	3.8
Maximum capacity (kg)	136





WA - 7

Walky®

Product features

- Foldable, Aluminum frame (easy to transport).
- Soft seat covered with eco-leather.
- Brakes with clamps.
- Four maintenance free 20 cm diameter wheels.
- Swivel front wheel.
- Adjustable height. • Ergonomic handle.
- Is pushed by third person. • Shopping basket,
- Accessibility of changing the direction while rollator.
- Possibility of using the footrest.

**Use:** Walking aid which improves stability and helps to transport the luggage especially useful when an injury or disability affects one leg only.



width (cm)	60	seat width	35,5
Length (cm)	70	Weight (kg)	8,2
Height (cm)	77-88	Maximum capacity (kg)	110
Length of legrest (cm)			90

WA - 8

Product features

- Steel frame.
- Soft armrest with easy adjustable of height.
- Easy adjustment of handle angle and height.
- Rear wheel with brakes.
- Padded platform helps people to move easily.

**Use:** Walking aid which improves stability and helps to maintain upright ambulation.



width (cm)	63
Length (cm)	73
Height (cm)	110-135
Maximum capacity (kg)	136



**Bubbly**<sup>®</sup>

**Bedy**<sup>®</sup>

**Electric Bed &  
Buble Mattress**



## PBM7

# Bubbly®

### Product features

- Model: PBM7.
- 130 Bubble pad design.
- Durable, medical grade PVC.
- Max patient weight 135kg ( 298lb ).
- Dimension: 200 x 90 x 7 cm / 79" x 35" x 3".
- Mattress weight 2.3kg ( 5lb ).
- Pump weight 1.4kg ( 3lb ).
- Operating power 220V - 50/60Hz.
- Provide superior treatment.
- Advance comfort therapy.
- Quiet, low vibration compressor.
- Easily attaches to bed frame.
- Adjustable pressure setting.



Adjustable hangers  
to suit any kind of  
hospital bed frame



High quality rubber tube

Free repairing kit



The most quiet pump

## P3FEB

### Product features

- 4-Functions electric hospital bed.
- Model: P3FEB.
- Size: 2125 x 970 x 430-830 mm.
- Max load 180kg ( 400lb ) in movement.
- Steel electrostatic spraying.
- Backrest tilting from  $0^{\circ}$  to  $75^{\circ}$  / (  $+,- 5^{\circ}$  ).
- Footrest tilting from  $0^{\circ}$  to  $40^{\circ}$  / (  $+,- 5^{\circ}$  ).
- Sit position function.
- Handset.
- Accessories: 4 Pieces 5" Silent Castors all with breakes, 1 pair Aluminum Side Rail, 2 Pieces ABS head and foot board, 4 Pieces Urine Hooks, 1 set MOTECK Motors, 1 Piece Electric Cable Holder, 4 Pieces Dust Cover.
- Option: I.V pole, Mattress, Over Bed Table, Back-up Battery, Beside Locker.



Back Rest Tilting  $\geq 70^{\circ}$



Hi-Low



Footrest Tilting



Sit Position



# Notes

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

# Partner<sup>®</sup>

*Let Us Improve Your Health*



[info@partner-medical.com](mailto:info@partner-medical.com)



[www.partner-medical.com](http://www.partner-medical.com)



Tel&Fax : +90 212 671 19 17



Ikitelli OSB MH. 6. CD.

Beyaz tower No:1 Office No:18

Basaksehir-Istanbul