

P - 100

Chairy®








Product features

- Frame Steel, single cross brace.
- Surface treatment powder coating.
- Front wheel 8"x1" front castor.
- 24" 9 spoke MAG wheel with grey PU tire.
- Padded flip back armrest with desk length.
- Swing away detachable footrest.
- Front fork Steel.
- Height adjustable Foot plate Plastic.
- Nylon seat and backrest.
- Two leg straps for safety.
- Foldable wheelchair.



Optional:

- Seat Belt - Elevating Legrests - Pneumatic Rear Wheels - Seat Cushion - Different Colors / Sizes

Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
						
40cm - 16" 46cm - 18" 50cm - 20"	40	92	66	112	19.8/17.8	115

Note : All dimensin in (cm , kg)