








P - 100 C

Chairy®

Product features

- Steel pediatric wheelchair.
- Yellow powder coating wheelchair frame.
- 20" spoke wheel with pneumatic tires.
- 6" front castor with PU tires.
- Powder coating steel hand rim.
- Padded flip back up armrest with desk length.
- Swing away detachable blue footrest with leg strip.
- Black powder coating aluminum alloy front fork.
- Height adjustable Plastic footplate.
- Nylon seat and backrest with sponge.
- Against backward overturn device.
- Foldable wheelchair.



Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 35cm - 14"	 35	 94	 59	 91	 15/20	 75

Note : All dimensin in (cm , kg)