

P - 103

Chairy®








Product features

- Powder coated grey Aluminum frame with single cross braces.
- 24" (61cm) 9-spoke MAG wheel with grey PU tire.
- 8x2" (20cm x 5cm) grey front castor with PU tires.
- Plastic front fork.
- Black plastic hand rim.
- Slope armrest with soft PU arm pad.
- Nylon seat with sponge.
- Swing in and away footrest.
- Height adjustable plastic footplate.
- Foldable wheelchair.



Optional:

- Seat Belt
- Pneumatic Rear Wheels
- Seat Cushion
- Different Colors / Sizes

Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 40cm - 16" 46cm - 18" 50cm - 20"	 46	 91	 61	 103	 21.5/19.0	 120

Note : All dimensin in (cm , kg)