








P - 105CP

Product features

- Ultra-light steel frame.
- Powder coating wheelchair frame.
- Easy to adjust chest holders to fix the user body upright position.
- Reclining high back and adjustable headrest.
- Stylish mag style wheels light up this product.
- The frame allows tilting between 3°-22° up to user demand.
- 16" spoke rear wheels.
- 6"x2" front castors.
- Powder coating steel front fork.
- Detachable armrest.
- Elevating footrest.
- Foldable wheelchair.



Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
 46 cm - 18 "	 40	 103	 61	 103	 22.2 / 20.4	 120

Note : All dimensin in (cm , kg)