


P - 340W

ENJOY THE INDEPENDENCE OF A POWER WHEELCHAIR WITH THE PORTABILITY OF A LIGHTWEIGHT MANUAL WHEELCHAIR.

Product features

Frame	Steel, Powder Coating
Controller	MAIKONG 50A
Motor	340W x 2 LINIX
Battery	Led - Acid 2x12V 40Ah
Charger	HP 24V, 4A
Max. speed	7 Km/h
Travel milage	33 Km
Footrest	Swing Away , Heel strap, Height Adjustable
Seat	Detachable Cusion Multi Size
Armrest	Flip Back Slope Armrest
Gradeability	10 Degree
Backrest	Half Folding Backrest
Front wheel	8" PU Castor
Rear wheel	12.5" Tire

Chairy®


Seat width


43cm - 17"
46cm - 18"
48cm - 19"

Seat depth


43

Wheelchair height


97

Wheelchair width


63

Wheelchair length


108

Weight
GW / NW


60/56

Maximum capacity


150

Note : All dimensin in (cm , kg)