

P - 250w

Chairy®

Enjoy the ultimate and durable driving experience with the P-250w most advanced with featured MAIKONG Control technology and innovative controls, easy to fold wheelchair.

Product features

Frame	Steel, powder coating
Controller	WUYANG 45A
Motor	250W x 2
Battery	2x12, 20 Ah
Charger	AC220 50~60Hz
Max. speed	6 Km/h
Travel milage	15 Km
Gradeability	6 degree
Armrest	Flip back slope armrest
Footrest	Swing in and away
Seat width	17" / 18" / 19"
Front wheel	10"
Rear wheel	18"
Side panel	Plastic
Other	With anti-tip



Seat width



43cm - 17"
46cm - 18"
48cm - 19"

Seat depth



43

Wheelchair height



97

Wheelchair width



45

Wheelchair length



78

Weight

GW / NW



60/56

Maximum capacity



120

Note : All dimensin in (cm , kg)