

CO - 2

CO - 3



Product features

- Foldable Chromed steel frame.
- Plastic commode seat.
- Fixed armrest and elevating footrest.
- Solid castors and rear wheels.
- Reclining high back with handle under the armrest.

Product features

- Foldable Chromed Steel frame.
- Plastic commode seat.
- Flip-up armrest and detachable footrest.
- Solid castors and solid rear wheels.

Note : All dimensin in (cm , kg)

	Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight GW / NW	Maximum capacity
CO-2	44 cm - 13"	43	129	64	117	28.7 / 25.2	100
CO-3	45 cm - 18"	43	90	60	106	21.4 / 19.6	100