

CO - 4



Product features

- Steel frame (Chrome paltded).
- leather upholstery.
- Flip down armrest.
- Detachable footrest.
- Paking brake mechanism.
- Plastic commode seat.

CO - 5



Product features

- Aluminum Frame.
- PE seat and back.
- round bucket.
- 3" wheels with parking mechanism.

Note : All dimensin in (cm , kg)

	Seat width	Seat depth	Wheelchair height	Wheelchair width	Wheelchair length	Weight	Maximum capacity
CO-4	44 cm - 17"	43	92	55	82	8.8	136
CO-5	44 cm - 17"	45	89	55	84	5	136