

PBM7

Bubbly®

Product features

- Model: PBM7.
- 130 Bubble pad design.
- Durable, medical grade PVC.
- Max patient weight 135kg (298lb).
- Dimension: 200 x 90 x 7 cm / 79" x 35" x 3".
- Mattress weight 2.3kg (5lb).
- Pump weight 1.4kg (3lb).
- Operating power 220V - 50/60Hz.
- Provide superior treatment.
- Advance comfort therapy.
- Quiet, low vibration compressor.
- Easily attaches to bed frame.
- Adjustable pressure setting.



Adjustable hangers
to suit any kind of
hospital bed frame



High quality rubber tube

Free repairing kit



The most quiet pump