

WA - 5

Walky®

Product features

- Made of Aluminum.
- Height adjustable. • Packed in pairs
- Armpit and hand parts covered in polyurethane foam, warm in contact.

Use:

Walking aid which improves stability and maintains upright ambulation.

Size	S	M	L
Height (cm)	84-114	114-134	134-154
Weight (kg)	1.6	1.7	1.8
Maximum capacity (kg)	115	115	115



WA - 6

Product features

- 3 Functions walker can be used as reciprocal, fixed or wheeled walker.
- Foldable light Aluminum frame.
- Stable and durable construction.
- Removable rubber walker tips.
- Additional removable front tips with 12 cm diameter.
- Ergonomic handles. • Height adjustable.

Use:

Walking aid which improves stability and helps maintains upright ambulation.

Height (cm)	77-96
Weight (kg)	3.8
Maximum capacity (kg)	136

