

CLAVICULA BANDAGE

Made of sponge that is covered with the terry fabric from the inside. Stabilize the clavicle thanks to the straps that pulls the shoulder to the back to give the appropriate position.


USAGE
Clavicle dislocations
Clavicle fractures



SU-2009

Size:

S 70-80 cm
M 80-90 cm
L 90-100 cm
XL 100-110 cm
2XL 110-120cm



SHOULDER SUPPORT (Neoprene)



Made of perforated neoprene fabric in 3 mm thickness. Provides warmth to the shoulder with protection in cases of the minor pain.

USAGE
Shoulder dislocations
After surgeries
Warming the Shoulder

SU-2010

Size:

S 34-38 cm
M 38-44 cm
L 44-50cm
XL 50-58 cm
2XL 58-66cm



SHOULDER SUPPORT (Wool)

Made of elastic cotton fabric specially designed to reduce sweating. Its inside is covered by smooth non-allergenic fabric. Designed to fit under clothes.

USAGE
Warming the Shoulder



SU-2011

Size:

S 34-38 cm
M 38-44 cm
L 44-50cm
XL 50-58 cm
2XL 58-66cm



SHOULDER ORTHOSIS



Made of thermoplastic covered by terry fabric from the inside. Has an optional elastic ligament. Stabilize the shoulder and the upper side of the arm.

USAGE
Stabilizing the shoulder and the upper part of the arm.
After surgeries.
Protection and Stabilization.

SU-2012

Size:

S 32-35 cm
M 36-39 cm
L 40-43 cm

