

PADDED ELBOW SUPPORT

Made of neoprene fabric that provides warmth and protection.
Used by athletes in cases of minor pains

USAGE

Protects elbow area against impacts.



SU-3003

Size:

S	22-25 cm
M	25-28 cm
L	28-31 cm
XL	31-34 cm
2XL	34-37cm



TENNIS ELBOW SUPPORT UNISIZE



Made of thermoplastic covered by terry fabric from the inside.

The angle can be adjusted easily
The Upper part can be adjusted in (0-15-30-45-60-75-90-105) degrees
The Lower part can be adjusted in (0-15-30-45-60-75-90-105) degrees

USAGE

Enables controlled mobilization of elbow.
medial and lateral epicondyle fractures.
After Surgeries.
elbow dislocation, elbow arthrolysis.

SU-3006

Size:

UNISIZE

ELBOW CONTRACTURE ORTHOSIS

Made of thermoplastic covered by terry fabric from the inside.

The angle can be adjusted easily
The Upper part can be adjusted in (0-15-30-45-60-75-90-105) degrees
The Lower part can be adjusted in (0-15-30-45-60-75-90-105) degrees

USAGE

Enables controlled mobilization of elbow.
medial and lateral epicondyle fractures.
After Surgeries.
elbow dislocation, elbow arthrolysis.



SU-3007

2 Size (S-M)

2 Side (Right-Left)

ELBOW IMMOBILIZER



Made of sponge that is covered with the terry fabric from the inside.

Stabilize the Elbow thanks to the metal pieces in it

USAGE

Protections at impacts.
Fractures.
Protection during training and sport activities.

SU-3008

Size:

UNISIZE