

## WRIST SPLINT LONG

Made of elastic cotton fabric specially designed to constrain the movements of the wrist thanks to the aluminum rod inside.

Provides better stabilization thanks to the elastic bandage on the wrist.

### USAGE

Provides immobilization and constrains the movements of the wrist.

Arthritis.

Rheumatism.



### SU-4001

#### Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25 cm



## WRIST SPLINT (NEOPRENE) LONG



Made of neoprene fabric that provides warmth and protection.

Constrains movements of the wrist thanks to the aluminum support inside.

### USAGE

Wrist tunnel cases.

Tendonitis, tenosynovitis.

After taking off gypsum.

After Surgeries.

### SU-4002

#### Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25 cm



## WRIST SPLINT (NEOPRENE) Right and Left

Made of neoprene fabric that provides warmth and protection.

Constrains movements of the wrist thanks to the aluminum support inside.

### USAGE

Wrist tunnel cases.

Tendonitis, tenosynovitis.

After taking off gypsum.

After Surgeries.



### SU-4004

#### Size:

UNISIZE

## WRIST SPLINT WITH THUMB SUPPORT LONG



Made of elastic cotton fabric specially designed to constrain the movements of the wrist thanks to the aluminum rod inside.

The thumb supporting rod stabilize it in the rest position.

Provides better stabilization thanks to the elastic bandage on the wrist.

### USAGE

Protection and stabilization of the wrist and the thumb.

Arthritis and rheumatism.

Wrist tunnel cases.

Pains in the thumb and wrist

### SU-4005

#### Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
XXL	23-25 cm

