

ADJUSTABLE WRIST SPLINT (NEOPRENE)

It is manufactured from perforated fabric. It is a type of splint that provides immobilization of the wrist in neutral position. It has an aluminium support inside and 2 plastic supports at the top. There is a bandage that provides stabilization in the laced system around the wrist and can be used in both hands.

USAGE

Carpal Tunnel Syndrome
In tendon and soft tissue injuries
Rehabilitation after surgical intervention



SU-4011

Size:
UNISIZE

NEOPRENE WRIST BANDAGE



It is a simple wrist bandage. It provides protection and heat to the wrist in case of wrist injuries, sprains and trauma. It improves blood circulation by maintaining the body temperature and thus speeds up the healing process.

USAGE

In wrist strains
To prevent sports injuries
In carpal dislocation

SU-4013

Size:
UNISIZE

NEOPRENE HAND WRIST SUPPORT

It is a simple wrist bandage that put thumb to inside. It provides protection and heat to the wrist in case of wrist injuries, sprains and trauma. It improves blood circulation by maintaining the body temperature and thus speeds up the healing process.

USAGE

Twists, Trumas, bruises.
Protection while doing sports.
Protection at work when the jobs involves using the wrist a lot.



SU-4014

Size:
UNISIZE

THUMB SUPPORT (DE QUEUVAIN)



It is a type of splint that immobilize thumb. The plastic atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

USAGE

CMC and MCP joint irritation
In de quervain syndrome
In ligament damage

SU-4015

Size:
UNISIZE