

THUMB SUPPORT

Made of thermoplastic material, It is a type of splint that immobilize thumb. The plastic atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

USAGE

CMC and MCP joint irritation
In de quervain syndrome
In ligament damage



SU-4016

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm



WRIST BAND WITH THUMB SUPPORT



It is a type of splint that immobilize thumb. The aluminium atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

USAGE

CMC and MCP joint irritation
In de quervain syndrome
In ligament damage

SU-4017

Size:

UNISIZE

WRIST BAND WITH THUMB SUPPORT LONG

It is a type of splint that immobilize thumb. The aluminium atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

USAGE

CMC and MCP joint irritation
In de quervain syndrome
In ligament damage



SU-4018

Size:

UNISIZE

WRIST SPLINT WITH DOUBLE SIDE



There is an anatomically shaped aluminium support that provides immobilization of the wrist in the neutral position. There is a bandage that increases the stabilization around the wrist. It is a standard size and can be personally adjusted with velcro.

USAGE

Carpal Tunnel Syndrome
In tendon and soft tissue injuries
Rehabilitation after surgical intervention

SU-4019

Size:

UNISIZE