

## WRIST SPLINT WITH THUMB SUPPORT NEOPRENE

There is an aluminium support under the wrist and near the thumb. It can use to both hands by changing aluminium support under the wrist. There is a velcro that increases the stabilization around the wrist.

### USAGE

In de quervain syndrome  
Carpal Tunnel Syndrome  
Rehabilitation after surgical intervention  
In wrist and thumb calcification  
In tendon and soft tissue injuries



**SU-4020**

**Size:**  
**UNISIZE**

## THERMOPLASTIC RELIEVING SPLINT



Hand and wrist orthosis manufactured from thermoplastic material. Keeps the fingers and the wrist in recovery position. Enables immobilization of wrist, used for resting position following plaster treatment of fractures and painful positions of the wrist.

### USAGE

For bone fractures of the wrist after plaster casting  
For resting in soft tissue injuries  
Carpal Tunnel Syndrome  
Rehabilitation after surgical intervention  
In tendon and soft tissue injuries

**SU-4021**

**Size:**

S	15-17 cm
M	17-19 cm
L	19-21 cm



## THERMOPLASTIC ANTI-SPASTICITY HAND SPLINT

Made of thermoplastic covered by terry fabric from the inside. Stabilizes the wrist and the fingers in a relaxing opened position.

### USAGE

Spasticity  
Wrist problems.



**SU-4022**

**Size:**

S	15-17 cm
M	17-19 cm
L	19-21 cm



## THERMOPLASTIC ANTI-SPASTICITY HAND SPLINT LUX



Made of thermoplastic covered by terry fabric from the inside. Stabilizes the wrist and the fingers in a relaxing opened position.

### USAGE

Spasticity  
Wrist problems.

**SU-4023**

**Size:**

S	15-17 cm
M	17-19 cm
L	19-21 cm

