

## THERMOPLASTIC DYNAMIC HAND-FINGER SPLINT

Made of thermoplastic covered by terry fabric from the inside.  
Pulls up the fingers and let them move a bit to help training them.

### USAGE

Rehabilitation.  
Extensor tendon.s.



**SU-4024**

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm



## THERMOPLASTIC KLEINERT HAND-FINGER SPLINT



Made of thermoplastic covered by terry fabric from the inside.  
Pulls up the fingers and let them move a bit to help training them.

### USAGE

Rehabilitation.  
Flexure tendon.

**SU-4025**

Size:

S	15-17 cm
M	17-19 cm



## MALLET FINGER SPLINT

Made of high quality plastic.  
Stabilize the upper part of the finger and prevent it from moving.

### USAGE

Mallet finger deformities.  
Interruption of tendons.  
Fractures and bruises.



**SU-5001**

Size:

S	4.5 - 5 cm
M	5 - 6 cm
L	6 - 7 cm
XL	7 - 8 cm
2XL	8 - 9 cm

## FROG SPLINT



Made of the Soft aluminum witch has the ability to reshape as wished.  
The inner side is covered with the plastoz witch helps stabilizing the finger.

### USAGE

Finger fractures.  
Strings and nerve injuries.  
Finger injuries and pain.

**SU-5002**

Size:

S	4 - 6 cm
M	6 - 8 cm
L	8 - 10 cm