

## POSTUREX BANDAGE

Made of elastic cotton fabric specially designed to reduce the Sweating.  
Has a plastic supporting reduces the ability to bend.  
Has side laces helps pulling the shoulders to keep the body's stabilization.  
It can be easily worn under the clothes.

### USAGE


Shoulder and back muscles pain.  
Other muscle problems.  
Long computer and car driving hours.  
Back bags carrying



**SU-6005**

### Size:

|     |            |
|-----|------------|
| S   | 70-80 cm   |
| M   | 80-90 cm   |
| L   | 90-100 cm  |
| XL  | 100-110 cm |
| XXL | 110-120cm  |



## MAGNETIC BACK CORSET



Adjustable straps design to fit all sizes the flexibility that gives you the freedom of movement. It has a magnetic belt that supports the spine to provide stability for those suffering from lower back pain


### USAGE

Great for lower back pain and supporting the back. It also keeps lower back in an upright or extended position it reduces pain in muscles, ligaments, discs and nerve roots of the lower back

**SU-6006**

### Size:

|     |            |
|-----|------------|
| S   | 70-80 cm   |
| M   | 80-90 cm   |
| L   | 90-100 cm  |
| XL  | 100-110 cm |
| XXL | 110-120cm  |



## ABDOMINAL CORSET

Made of elastic cotton fabric specially designed to reduce sweating.  
Stabilize the abdomen.  
Can be adjusted as needed.  
Reduces the the pains in the abdomen area.  
Designed to fit under clothes.

### USAGE


Post abdominal surgery and birth  
Peritoneal injury  
Abdominal sagging  
Post umbilical and ventral hernia surgery



**SU-6007**

### Size:

|     |            |
|-----|------------|
| S   | 70-80 cm   |
| M   | 80-90 cm   |
| L   | 90-100 cm  |
| XL  | 100-110 cm |
| XXL | 110-120cm  |



## CHEST CORSET



Made of elastic cotton fabric specially designed to reduce sweating.  
Its inside is covered by smooth non-allergenic fabric.  
Designed to fit under clothes.

### USAGE

Sternal Support after coronary surgery  
Costa and sternum fractures  
Post chest surgery and trauma

**SU-6008**

### Size:

|     |            |
|-----|------------|
| S   | 80-90 cm   |
| M   | 90-100 cm  |
| L   | 100-110 cm |
| XL  | 110-120 cm |
| XXL | 120-130cm  |

