

(NEOPRENE) KNEE SUPPORT CLOSED PATELLA

Made of neoprene fabric
 Contains pores that help for ventilation.
 Support knee joint, reduce pain by relieving the pressure on the joint of the knee.
 Sustains body's own release of heat effectively improve blood circulation, promote metabolism, enhance immunity also applies compression, prevents edema.

USAGE

Arthritis. Intra articular liquid increase
 Pain induced by knee cap
 Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities
 When heat and compression is desired



SU-7001

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



(NEOPRENE) KNEE SUPPORT CLOSED WITH SILICONE PAD



Made of neoprene fabric
 Contains pores that help for ventilation.
 Support knee joint, reduce pain by relieving the pressure on the joint of the knee.
 The perimeter of the patella is closed and supported with a silicone pad.
 Sustains body's own release of heat effectively improve blood circulation, promote metabolism, enhance immunity also applies compression, prevents edema.

USAGE

Arthritis. Intra articular liquid increase
 Pain induced by knee cap
 Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities
 When heat and compression is desired

SU-7002

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



(NEOPRENE) KNEE SUPPORT PATELLA SUPPORT

Made of neoprene fabric
 Contains pores that help for ventilation.
 Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Can be used daily as it does not compress your knee cap.
 Recommended for people who have pain in their knee cap with and difficulty when climbing stairs or squatting.

USAGE

Arthritis. Intra articular liquid increase
 Pain induced by knee cap
 Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities
 When heat and compression is desired



SU-7003

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



(NEOPRENE) KNEE SUPPORT PATELLA WITH LIGAMENT SUPPORT- SHORT



Made of neoprene fabric
 Contains pores that help for ventilation.
 Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Flexible wire on the sides increase stabilization and decrease the pressure on the knee.

USAGE

Arthritis, bursitis, tendinitis
 Sprains and stains
 Meniscus problems
 Post-op or trauma for protection purposes
 Ligament instability
 Sporting injuries

SU-7004

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

