

KNEE SUPPORT WITH HINGED STABILIZING JOINTS- FULL OPEN FRONT

Made of neoprene fabric
Contains pores that help for ventilation.
Aluminum joints on both sides provide maximum support, allow controlling the motion and prevent knee hyperextension.
Velcro fasteners on the upper and lower calves increase stabilization.
Velcro fasteners on the upper and lower calves increase stabilization.
Since it can be completely open from the front side, provides ease of use to the elderly, heavy weight and post-op patients.

USAGE

Ligament injury and instability
Meniscus lesions
1 ve 2 degree sprains and stains
Post-op or trauma for protection purposes



SU-7009

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



KNEE SUPPORT PATELLA WITH CROSSED LIGAMENT SUPPORT-NEOPRENE



Made of neoprene fabric
Contains pores that help for ventilation.
Contains opened the patella supported by a pad. Medial (inside) and lateral (outside) Flexible wire offer overall knee joint stability.
Transverse straps adhering with 3-point principle provide stabilization for ant. Cruciate ligaments and patella.

USAGE

Ant. Cruciate ligament injury and knee instability
Patellar stabilization
Meniscus lesions
Post-op or trauma for protection and rehabilitation purposes

SU-7010

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



KNEE SUPPORT FOR PATELLAR TENDON

Made of neoprene fabric that provides warmth and protection.
The special pad supports the patellar tendon and with the applied compression reduces pain.

USAGE

Os-good Schlatler's disease
Patellar tendinitis
Patellar chondromalacia
Patellofemoral pain



SU-7011

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



PATELLAR TENDON BAND UNISIZE



Made of neoprene fabric that provides warmth and protection.
The special pad supports the patellar tendon and with the applied compression reduces pain.
Does not restrict joint motion.
Providing comfortable all-day wear.
The knee support strap is made with neoprene so it provides a comfortable and flexible fit.
Used for both knees.

USAGE

Os-good Schlatler's disease
Patellar tendinitis
Patellar chondromalacia
Patellofemoral pain

SU-7012

**Size:
Standard**