

THIGH SUPPORT

Made of neoprene fabric
Contains pores that help for ventilation.
It supports thigh's muscles
effectively improve blood circulation,
Keeps the heat, reduces the pressure,
reduces the severity of muscle's pain

USAGE

Muscle injuries to the thigh
Thigh pain
Protection during and after the sportive
activities
When heat and compression is needed



SU-7027

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



CALF SUPPORT



Made of néoprène fabric
Contains pores that help for ventilation.
It supports thigh's muscles
effectively improve blood circulation,
Keeps the heat, reduces the pressure, re-
duces the severity of muscle's pain

USAGE

Muscle injuries to the calf
Calf pain
Protection during and after the sportive
activities
When heat and compression is needed

SU-7028

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



HIP SUPPORT

It is manufactured to support thigh. It re-
duces pain by radiate heat to thigh area.
Velcros on it provides easy usage. Include
belt and so it provides comfortable usage
to patients for daily usage.

USAGE

Muscle injuries to the thigh
Thigh pain
Protection during and after the sportive ac-
tivities
When heat and compression is needed



SU-7029

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm



HIP ABDUCTION ORTHESIS



Hip abduction Orthosis anatomically de-
signed and manufactured of thermoplastic
material, providing hip stabilization.
Parts in contact with the skin are coated
with terry cloth. Flexion adjustment possi-
ble at every 15 degrees, starting from 15
up to 120 degrees.

Abduction angle is adjustable.

USAGE

Stabilization of total hip prosthesis or revi-
sion
Femoral head resection
Femoral or acetabular osteotomy
Recurrent hip dislocation

SU-7030

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47cm

