

ROM WALKER

It comes with 3 pads (12.5°) so it can be adjusted as needed.
Physiological walk due to special sole.
Plantar flexion: 0°, 10°, 20°, 30°, 40°
Dorsiflexion: 20°, 10°, 0°

USAGE

After ligament, soft tissue and tendon surgeries.
For post operative and follow up treatments of the talus, calcaneus and medial malleolus fractures.
For early functional therapy after conservative or surgical treatment for rupture for the Achilles tendon



SU-8009

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



ACHILL BOOT



It comes with 3 pads (12.5°) so it can be adjusted as needed.
Physiological walk due to special sole.

USAGE

After ligament, soft tissue and tendon surgeries.
For post operative and follow up treatments of the talus, calcaneus and medial malleolus fractures.
For early functional therapy after conservative or surgical treatment for rupture for the Achilles tendon

SU-8010

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



REFLEX AFO

Ankle orthosis made of plastic material.
Provides static dorsiflexion assistance and lateral stability for entire ankle/foot area.
Provides ease of use with velcro fasteners going around the ankle.
Use a scissor to trim insole for exact shoe fit.
May be modified with a heat gun.

USAGE

Foot-drop deformity caused by nerve injuries such as hemiplegia, lumbar disc hernia
When 90 degrees ankle angle is needed



SU-8011

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



FOOT GUARD



T shaped ankle orthosis made of thermoplastic material with terry cloth coating.
Ensures the fix 90 degrees stand, prevent rotation of the ankle with T Shape.

USAGE

Stabilizing the position of the leg on hospital patients

SU-8012

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm

