

HALLUKS VALGUS NIGHT SPLINT

Foam laminated fabrics, the inner surface is made of terry cloth. Position of the thumb to the desired position is hidden inside the fabric with aluminum splint can be adjusted. designed to be used at Night (in sleep) .

USAGE

Helps to keep the big toes aligned, designed to use during night time.



SU-8013

Size:
Standard

HALLUKS VALGUS DAY SPLINT

Velcro and elastic bandage is made of neo-prene fabric. Day (shoes, slippers ... in) designed to be used at night.

USAGE

Helps to keep the big toes aligned, designed to use during day time.



SU-8014

Size:
Standard

CASTHOE

with adjustable Velcro tapes can also be applied to both feet.

USAGE

Used instead of shoes in case of plaster feet



SU-8015

Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm

