

# CERVICAL COLLARS

## CERVICAL COLLAR (Sponge)

Made of compressed Sponge witch is covered by cotton fabric.  
Keeps the neck warm and stabilize the neck and the head.

### USAGE

Muscles pain.  
Neck arthritis.  
Sliding discs.  
After the traumas and bruises



**SU-1001**

### Size:

S 34-36 cm  
M 36-38 cm  
L 38-40 cm  
XL 40-42 cm



## NELSON COLLAR



Made of the foam witch has holes in it that helps coaling the neck, allows the x-ray to goes throw it, non-allergic, reduces the pain by reducing the ability of the head and the neck to move.

### USAGE

Muscles pain.  
Neck arthritis.  
Sliding discs.  
After the traumas and bruises

**SU-1002**

### Size:

S 34-36 cm  
M 36-38 cm  
L 38-40 cm  
XL 40-42 cm



## PHILEDELPHIA COLLAR

Made of the foam, reduces the ability of moving the neck and head thanks to the plastic supports on the front or the back, allows the x-ray to go throw it.

### USAGE

Soft tissue trauma.  
Ligaments and bones injuries.  
Sliding discs.  
Stabilizations in the cases of the uncon-  
sciousness.



**SU-1003**

### Size:

S 34-36 cm  
M 36-38 cm  
L 38-40 cm



## PHILEDELPHIA COLLAR (Front Open)



Made of the foam, reduces the ability of moving the neck and head thanks to the plastic supports on the front or the back, allows the x-ray to go throw it, the front has a hole in the Trachea area to be used in case of emergency.

### USAGE

Soft tissue trauma.  
Ligaments and bones injuries.  
Sliding discs.  
Stabilizations in the cases of the uncon-  
sciousness.

**SU-1004**

### Size:

S 34-36 cm  
M 36-38 cm  
L 38-40 cm  
XL 40-42 cm

