

## ELBOW SUPPORT KIDS

Made of neoprene fabric that provides warmth and protection.  
Thanks to the three-headed silicone pad inside, it applies constant pressure and compression to the epycondil heads.

### USAGE

Lateral epicondylitis.  
Internal articular fluid.  
Elbow pain



**SU-9009**

### Size:

S	16-22 cm
M	22-28 cm
L	28-34 cm



## WRIST BANDAGE KIDS



manufactured of neoprene fabric  
Provides protection for the wrist, sprains and bruises wrist  
Retains body heat to improve bloody circulation and accelerate wound healing

### USAGE

Twists, Trumas, bruises.  
Protection while doing sports.  
Protection at work when the jobs involves using the wrist a lot.

**SU-9010**

**Size :**  
**Standard**

## WRIST SPLINT KIDS

Made of neoprene fabric that provides warmth and protection.  
Constraints movements of the wrist thanks to the aluminium support inside.

### USAGE

Wrist tunnel cases.  
Tendonitis, tenosynovitis.  
After taking off gypsum.  
After Surgeries.



**SU-9011**

### Size:

S	9-11cm
M	11-13 cm
L	13-15 cm



## WRIST SPLINT WITH THUMB SUPPORT (KIDS )



There is an aluminium support under the wrist and near the thumb. It can use to both hands by changing aliminium support under the wrist. There is a velcro that increases the stabilization around the wrist.

### USAGE

Wrist tunnel cases.  
Tendonitis, tenosynovitis.  
After taking off gypsum.  
After Surgeries.

**SU-9012**

### Size:

S	9-11cm
M	11-13 cm
L	13-15 cm

