

## FIGURE OF EIGHT BANDAGE

Made of neoprene fabric  
Contains pores that help for ventilation.  
The Crossed lacing provides extra pressure and stabilization thanks to the rubber pads on the sides.

### USAGE

Ligament instability and rapture  
Edema prevention  
Post-op or trauma for protection and rehabilitation purposes



**SU-9017**

**Size :**  
**Standard**

## PAVLIK SUPPORT



Made of 3 cm cotton belts, non-elastic, stabilize the hip and the feet in a perfect position.

### USAGE

Congenital Dislocation

**SU-9018**

**Size:**

<b>S</b>	<b>0-3 Month</b>
<b>M</b>	<b>3-6 Month</b>
<b>L</b>	<b>6-9 Month</b>

## FREJKA ORTHOSIS

Manufactured from 12 mm top quality molded plastazote foam.  
Holes allow air to circulate, soothing support and relief.

Latex-free, non-toxic and hypoallergenic plastizode foam material will not cause skin irritation. Velcro straps are used for easy donning and doffing. Plastizode foam material also is radiolucent which allows X-ray, it possible to use CT, MRI. Positions at the hip joint in flexion and abduction.

### USAGE

Hip subluxations



**SU-9019**

**Size:**

<b>S</b>	<b>0-3 Month</b>
<b>M</b>	<b>3-6 Month</b>
<b>L</b>	<b>6-9 Month</b>

## UMBILICAL HERNIA TRUSS



made of flexible cotton fabric. The pad applies pressure to the inguinal hernia. Perineal strap prevents sliding. It comes in standard size but adjustable.

### USAGE

UMBILICAL HERNIA

**SU-9020**

**Size :**  
**Standard**