

## SILICON BREAST

It is used to provide a natural appearance when chest tissue is removed after mastectomy, to reduce pain caused by muscle imbalances, to protect the surgical wound from external bumps and crushes. The laundry is made of cotton and non-allergenic special fabric and silicone material. Shoulder straps are softened for comfortable use.

The pockets on both sides allow the prosthesis to be used as right and left.



**SU-100**

**Size:**  
**Standard**

## SILICON BREAST'S BRA



The prosthetic bra, which is very comfortable and stylish to wear, is specially designed for breast prosthesis. It supports the prosthesis and weight and provides a sense of confidence and comfortable use. It is made of cotton fabric, the inner surface of which does not irritate the skin. Shoulder straps are softened for comfortable use. The pockets on both sides allow the prosthesis to be used both left and right.

**SU-101**

**Size:**  
**Standard**

## PORTABLE CLOSED SEAT LIFTER (WITHOUT COVER)

Made of 3 cm cotton belts, non-elastic, stabilize the hip and the foot in a perfect position.

**USAGE**  
Congenital Dislocation



**SU-102**

**Size:**  
**Standard**

## PORTABLE CLOSED SEAT LIFTER (WITH COVER)



Latex-free, non-toxic and hypoallergenic plastizode foam material will not cause skin irritation. Velcro straps are used for easy donning and doffing. Plastizode foam material also is radiolucent which allows X-ray, it possible to use CT, MRI. Positions at the hip joint in flexion and abduction.

**USAGE**  
Hip subluxations

**SU-103**

**Size:**  
**Standard**

## SILICON EXERCISE BALLS

Made of high quality silicone material which does not cause allergies and has no smell. Designed to withstand the pressure applied on, without losing their shape.

Available in three different colors pressure

**USAGE**

- Used so as to rehabilitation.
- Three different hardnesses existing.
- Red – Hard, Blue – Medium, Green – Soft



**SU-104**

**Size:**  
**Standard**