

STEEL INSOLE



made of stout leather and supported with stainless steel bridge has a longitudinal and metatarsal support.

Distributes the standing body weight evenly and reduces load and stress on the joints. Rigid surface support does not allow the longitudinal arc to collapse.

USAGE

Pes planus-planovalgus
Pes cavus
Transverse arc support

SU-8035

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



SILICONE METATARSAL PAD WITH RING

It is made of fabric material that does not harm the skin. The position of the silicone ring is adjusted by inserting it on the 2nd finger. Metatarsals relieve pain by removing the pressure on their heads.

USAGE

Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.



SU-8118

Size :
Standard

3/4 INSOLES (SILICONE)



Middle degree of hardness, is made from 100% silicon. 3/4 in size and offers comfortable use.

USAGE

Provides anatomical support to the sole and enables even distribution of pressure on the sole. Has a relaxing and comforting effect.

SU-8018

Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm



SILICON HEEL CUP

made of high quality of silicone material. Specially designed to fit figure of the foot. Softened calcaneal epin area lessens the load on the heels which causes pain. Lessens the pain at foot and calf.

USAGE

Heel spur
Plantar fasciitis
Avescular necrosis of calcaneus (sever's disease)
Heel pain



SU-8019

Size:

S	34-37 cm
M	38-39 cm
L	40-41 cm



SILICON TOE SEPERATOR



made of high quality of silicone material. Comfortable for daily use. Prevents deviation of the toe to the lateral in hallux valgus deformity.

USAGE

Helps proper positioning of toes in the shoes and prevents overlapping.

SU-8020

Size :
Standard

FOOT CARE