

## CERVICAL COLLAR WITH / WITHOUT CHIN SUPPORT (Plastic)

Made of PVC with chin support in a design that can be modified as wanted to stabilize the neck and reduce the pressure on the spine

**USAGE**  
Muscles pain.  
Neck arthritis.  
Sliding discs.  
Cases of radiculopathy  
After the traumas and bruises



**SU-1005 / SU-1006**

**Size:**  
S 34-36 cm  
M 36-38 cm  
L 38-40 cm  
XL 40-42 cm



## WIRE-FRAME COLLAR JR



Made of the aluminum covered by polyurethane to stabilize the neck as wished.

**USAGE**  
Discs hernia.  
Arthritis.  
Recovery from the surgery.  
Neck injuries.  
Tissue injuries

**SU-1007**

**Size:**  
S 34-36 cm  
M 36-38 cm  
L 38-40 cm  
XL 40-42 cm



## EMERGENCY NECK COLLAR

Made of the PVC plastic with a foam layer on the inner face of it.  
Can be adjusted to 3 different sizes.  
Has a lock system that helps stabilizing the neck.  
Has a hole in front of it to provide faster access to the neck in cases of emergency.  
Allows the X-ray to go through it.

**USAGE**  
Stabilizing the neck during the emergencies.  
In coma cases.



**SU-1008**

**Size:**  
S 34-36 cm  
M 36-38 cm



## MINERVA ORTHOSIS



Made of the Thermoplastic which is covered from the inside by sponge that helps stabilizing the neck and the upper part of the body.  
Stabilize the chest and reduces the neck's ability of movement.

**USAGE**  
After the surgeries of the muscles and ligaments.  
Arthritis.  
Bone fractures.

**SU-1009**

**Size:**  
S 34-36 cm  
M 36-38 cm  
L 38-40 cm

